# velocity®

*charge it up*

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
*This diagram is to help familiarize you with the unique parts of the Velocity. These parts are referenced in the instructions and used when speaking with a technical service representative.*

## KEY COMPONENTS

1. Integrated Pyramid
2. Housing
3. Posterior Toe Plate
4. Anterior Toe Plate
5. Posterior Toe
6. Anterior Toe
7. Foreheel
8. Heel Pad
9. Toe Pad
10. Heel Wedge (optional)
11. Toe Wedge (optional)

## PACKAGE CONTENTS

- (1) Velocity Foot
- (1) Foot Shell
- (1) CPI Sock
- (1) Velocity Wedge Kit

## TOOLS REQUIRED

- (1) 4mm Hex Key

## TOOLS RECOMMENDED

Foot Horn

## PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

## MOUNTING

Use only high quality proximal endoskeletal components.

## ASSEMBLY AND DISASSEMBLY (for sock replacement)

Use the FootHorn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

## STATIC ALIGNMENT

For optimal function have the patient's weight balanced between the forefoot and heel. The load line divides the foot at 1/3 heel lever.

## DYNAMIC ADJUSTMENTS

| DESIRED RESULT              | ALIGNMENT CHANGE                                     | COMPONENT CHANGE    |
|-----------------------------|------------------------------------------------------|---------------------|
| <b>Firmer Toe Response</b>  | Plantarflex the Velocity or move load line posterior | Install Toe Wedge*  |
| <b>Softer Toe Response</b>  | Dorsiflex the Velocity or move load line anterior    | Remove Toe Wedge*   |
| <b>Firmer Heel Response</b> | Dorsiflex the Velocity or move load line anterior    | Install Heel Wedge* |
| <b>Softer Heel Response</b> | Plantarflex the Velocity or move load line posterior | Remove Heel Wedge*  |

*\* To install or remove toe and heel wedges, refer to the Velocity Wedge Kit Instructions. Pre-installed toe wedges cannot be removed.*

## GAIT MATCHING® GUIDELINES

Contact College Park Technical Service if you have questions about modifying the gait match.

## WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

## **WARRANTY INSPECTION / MAINTENANCE INFORMATION**

College Park recommends that you schedule your patients for check-ups according to the warranty Inspection schedule below.

*High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.*

- Composites and Adapters
- CPI Sock
- Foot Shell

### **WARRANTY INSPECTION SCHEDULE FOR VELOCITY: ANNUALLY.**

## **TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365**

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST).


After hours, an emergency Technical Service number is available to contact a College Park representative.

## **LIABILITY**

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## **CAUTION**

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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©2019 College Park Industries, Inc. All rights reserved.