friction brake knee

ركبة ممانعة للاحتكاك • knæprotese med friktionsbremse • Reibungsbremsknie
 γόνατο με πέδη τριβής • rodilla con freno de fricción • kitkajarrupolvi
 genou avec frein à friction • ברך בלם חיכוך • ginocchio con freno ad attrito
 knie met wrijvingsrem • kne med friksjonsbrems • kolano z hamulcem ciernym
 joelho com travão de fricção • joelho com freio de fricção
 коленный сустав с фрикционным тормозом • kolenný kĺb s trecou brzdou
 friktionsbromsande knä • sürtünme frenli diz • 摩擦制动膝盖

college park
 TECHNOLOGY for the HUMAN RACE

technical instructions

التعليمات الفنية • Teknisk vejledning • Technische Anleitung • Τεχνικές οδηγίες • Instrucciones técnicas
 Tekniset ohjeet • Instructions techniques • הוראות טכניות • Istruzioni tecniche • Technische instructies
 Tekniske instruksjoner • Instrukcje techniczne • Instruções técnicas • Instruções técnicas
 Технические инструкции • Technické pokyny • Tekniska anvisningar • Teknik Talimatlar • 技术说明

FIGURE 1


FIGURE 2


FIGURE 3


FIGURE 4



FIGURE 5



FIGURE 6




FIGURE 7


FIGURE 8



FIGURE 9

PACKAGE CONTENTS

- (1) College Park Guardian Knee (1) Knee External Lock Kit

This diagram is to help familiarize you with the unique parts of the Guardian Knee. These parts are referenced in the instructions and used when speaking with a technical service representative.

KEY COMPONENTS (FIGURE 1)

- | | | |
|-------------------------------|---|-----------------------------|
| A. Pyramid | B. Upper Assembly | C. Lower Assembly |
| D. 30mm Pylon Receiver | E. 4mm Clamp Screw (Torque 10 N-m) | F. Release Lever |
| G. Toe Load Screw | H. Lock Release Screw | I. External Lock Kit |

PRODUCT DESCRIPTION

The Guardian Knee is constructed with an integrated pyramid (proximal) and 30 mm pylon receiver (distal).

INTENDED USE

The Guardian Knee, intended for transfemoral amputees, is a prosthetic device designed to restore some function of an anatomical knee joint.

INDICATIONS:

Lower limb amputations above the knee

CONTRAINDICATIONS:

None known

TECHNICAL SPECIFICATIONS

Material	Aluminum	Assembly Weight	533 g
Joint Type	Single-Axis	Patient Weight Limit	275 lb (125 kg)
Flexion	145°	Warranty	2 years
Build Height	1.0 in (2.6 cm)	Torque (Clamp Screw)	10 N-m

BUILD HEIGHT (FIGURE 2)

A	Overall height	5.2 in (13.1 cm)
B	Dome to knee center	1.0 in (2.6 cm)
C	Dome to tube end contact	2.6 in (6.6 cm)

BENCH ALIGNMENT (FIGURE 3)

- A.** Determine the heel height **B.** Determine the socket flexion **C.** The load line bisects through the pylon

STATIC ALIGNMENT

With the alignment reference line through the pylon, plantarflex or dorsiflex the foot until the load line is balanced between 1/3 heel and 2/3 toe lever.

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DYNAMIC ADJUSTMENTS

All adjustments can be made using a 4mm Allen wrench.

Note: Patient must be sitting when adjustments are made.

Note: The effect of a dynamic adjustment can be felt with as little as 1/8 - 1/4 turn.

Caution: Flexion and extension must be possible at all settings.

STANCE PHASE

LOAD DEPENDENCE (BRAKE) (FIGURE 4)

Factory setting = Screw is at minimum (most brake sensitivity).

Controls the brake sensitivity. When load is applied, the knee will not bend until the load is displaced. Turning the load adjustment clockwise will increase the amount of load required to initiate braking.

Note: The factory setting for load dependence should generally not require adjustment. However, it may be necessary in cases such as patients with higher body weights.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Brake function too strong; Swing too difficult to initiate	Increase Load Dependence (brake less sensitive)	Turn L Clockwise	
Brake function not sufficient; Knee too unstable during stance	Decrease Load Dependence (brake more sensitive)	Turn L Counterclockwise	

STANCE FLEXION ADJUSTMENT (FIGURE 5)

Factory setting = Turn screw clockwise to the point where resistance is felt, then reverse about 1/8 turn (counterclockwise).

Controls the amount of stance flexion. Adjustments can be made to alter the amount of stance flexion before the brake is activated.

Note: Stance flexion adjustment works together with load dependence to control the brake function. If decreasing the stance flexion adjustment, it may be necessary to also increase the load dependence adjustment.

Caution: Do not over tighten. May lead to brake sticking which can cause unstable gait.

Caution: Do not under tighten. Decreasing stance flexion may result in reduced brake engagement if load dependence is not also adjusted.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Not enough flexion	Increase Stance Flexion	Turn S counterclockwise	S ↗ +
Too much flexion	Decrease Stance Flexion	Turn S clockwise	- ↘ S

SWING PHASE

EXTENSION ASSIST (FIGURE 6)

Factory setting = Screw is at midpoint of range.

Controls the rate of extension assist.

Note: Knee must be fully flexed to access extension adjustment.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Extension too slow or Excessive heel rise	Increase Extension Assist	Turn E Clockwise	+ ↖ E
Extension too fast or Heel rise not sufficient	Decrease Extension Assist	Turn E Counterclockwise	E ↗ -

EXTERNAL LOCKING OPTION

Factory Setting = Stance Control Mode

The College Park Guardian can be converted to a locking knee joint, allowing the patient to control when the knee is locked or unlocked. To achieve this, the lock must be enabled, and the external lock kit must be installed. Refer to External Lock Kit Fabrication Instructions for detailed information and mounting instructions.

Stance Control Mode

Lock is disabled; Knee utilizes the brake mechanism during stance

Locking Mode

Lock is enabled; Allows for manual locking and unlocking of knee

ENABLING LOCKING MODE (FIGURE 7)

Note: Do not remove the release screw or apply Loctite.

1. Turn the release screw counterclockwise, about 3 ¼ turns.
2. Continue with installation of the External Lock Kit.

Caution: Test the knee to ensure it does not flex, prior to patient fitting.

DISABLING LOCKING MODE

Turn the release screw clockwise, about 3 ¼ turns. The knee will now be in stance control mode.

Caution: Test the knee to ensure it swings freely without weight applied, prior to patient fitting.

USING THE MANUAL LOCK FEATURE (FIGURE 8)

Caution: Test the dual function lever after installation and all adjustments to ensure proper function.

Note: Instruct the patient on proper use of the dual function lever. Make sure they understand the operating positions. Visit www.college-park.com to find fabrication instructions.

BENCH ALIGNMENT (FIGURE 3)

A. LOCKED POSITION

Standing/Walking: The knee will not flex.

B. LOCK RELEASE POSITION

Pull upward on the lever only half-way, without passing the detents. Releasing the lever will return it to the locked position.

C. UNLOCKED POSITION

Pull upward on the lever all the way, passing the detents. There will be a “click” sound and feeling. Push down on the lever to return to the locked position.

TO SIT DOWN

1. Patient should pull upward on the dual function lever to allow the knee to flex.
2. If toe load feature is enabled: to unlock the knee, the patient must shift their weight forward slightly (apply a toe load or hyperextend knee) while also pulling upward on the lever.


LOCK RELEASE SENSITIVITY (FIGURE 9)

The Guardian may be adjusted to require an extension moment (toe load) to unlock the knee. In this case, pulling on the dual function lever alone will not unlock the knee. If this feature is desired, perform the following:

1. Turn the TL screw clockwise, about 3 ½ turns.
2. To unlock the knee, the patient must shift their weight forward slightly (apply a toe load or hyperextend knee) while also pulling upward on the lever.

ADJUSTING LOCK RELEASE (TOE LOAD)

Factory setting = Screw is at minimum (NO toe load required).

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Unlocking too easy; requires too little toe load	Increased toe load response	Turn TL clockwise	
Unlocking too difficult; requires too much toe load	Decreased toe load response	Turn TL counterclockwise	

WARNINGS

- Flexion and extension must be possible at all settings.
- Avoid pinching hazards! Do not place fingers near the flexing area of the knee.
- Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to noise, sudden loss of function, brake release sticking, etc.
- Do not disassemble the knee. Contact College Park to arrange a repair or replacement.
- Do not expose this product to moisture, i.e. fresh water, corrosive materials, salt water or pH extremes. If the knee encounters moisture, wipe it dry using a lint-free cloth.
- Contaminants such as dirt and the use of lubricants or powder may affect the function of the knee brake and lead to failure.
- Do not use compressed air to clean the knee as it can push dirt inside of the knee.
- If the locking function is utilized, verify that the manual lever and nylon lanyard are not obstructed by the cosmesis.
- If modifying the prosthesis, ensure that the nylon lanyard is secure to prevent damage from grinding tools.

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Knee Assembly, External Lock

WARRANTY INSPECTION SCHEDULE FOR COLLEGE PARK GUARDIAN: SIX MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

*“Serious incident” is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

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©2022 College Park Industries, Inc. All rights reserved.