


ORDER FORM FOR ALLELES[™] LiL Covers

Order Date: Company Name:		Original PO#:	Billing PO#:
		Prosthetist:	
Patient Name:		Account Number:	
Email:		Telephone:	Fax:
1. SELECT DESIGN:		SHAPE: ase see 'SHAPE' diagrams	3. LEFT OR RIGHT:
LiL Alpine		following page.	Left Leg
LiL Baller LiL Nouveau LiL Stea		AK	Right Leg
LiL Deener LiL OM LiL Victor		BK1 BK2	
		ВКЗ	
Please reference the measuring diagrams on following page and be sure to select the unit of measure.			
4. COVER LENGTH (A): 5. TOP WIDTH (B): 6. BOTTOM WIDTH (C):			
inches inches inches			
millimeters millimeters millimeters			
7. MAIN COLOR		8. DETAIL COLOR	
14K Gold Cobalt Purple	Teal	14K Gold Coba	lt Purple Teal
Army Green Forest Green Racecar Yellow	White	Army Green Fores	st Green 🔲 Racecar Yellow 🗌 White
Baby Blue Grey Green Red	Skin #1		Green Red Skin #1
Baby Pink Lime Green Rose Gold Baby Yellow Maroon Seafoam Green	Skin #2	Baby Pink Lime	Green Rose Gold Skin #2 on Seafoam Green Skin #3
Black Navy Silver	Skin #4	Black Navy	
Charcoal Orange Spanish Gold	Skin #5	Charcoal Oran	ge Spanish Gold Skin #5
9. STRAP COLOR ANYTHING ELSE WE SHOULD KNOW?:			
Black (Default. No charge.)			
Tan (Additional Charge.)			
Red (Additional Charge.)			



SHAPE





Measure from the apex of the mechanical knee to the top of the footshell. A Leave a 1/2 inch (13mm) gap between the top of the footshell and bottom of the cover.



BK2

BK1 Length

Measure from the lowest dip in the front of the socket to the top of the footshell. 🗙 Leave a 1/2 inch (13mm) gap between the top of the footshell and bottom of the cover.

(B) AK Top Width


MEASUREMENTS

AK Bottom Width

(C)


Measure the calf circumference of the sound limb.

BK1 Top Width




 \bigstar Measure the circumference of the socket 1 inch (25mm) below the dip in the front of the socket.

BK2 Top Width



4 Measure the circumference of the socket 1 inch (25mm) below the dip in the front of the socket.

BK3 Top Width






BK1 Bottom Width




Measure the circumference right above the ankle bone of the sound limb.

BK2 Bottom Width



☆☆☆ Measure the circumference of the socket 2 inches (51mm) above the top of the footshell.

BK3 Bottom Width




Measure the widest circumference near the bottom of the socket. We will contact you for additional measurements



Measure from the lowest dip in the front of the socket to the top of the footshell. 🖈 Leave a 1/2 inch (13mm) gap between the top of the footshell and bottom of the cover.





Measure from the lowest dip in the front of the socket to the top of the footshell. \bigstar Leave a 1/2 inch (13mm) gap between the top of the footshell and bottom of the cover.

 \bigstar Measure the circumference of the socket 1 inch (25mm) below the dip in the front of the socket.

ALLELES LiL Covers



LIL ALPINE (SHOWN IN SKIN #5 WITH BABY BLUE DETAILS)



LIL BALLER (SHOWN IN BLACK WITH RED DETAILS)



LIL DEENER (Shown in Baby Pink with 14K Gold Details)



LIL GLIDE (SHOWN IN BABY YELLOW WITH RACECAR YELLOW DETAILS)



LIL HIFI (Shown in Baby Blue with Orange Details)



LIL NOUVEAU (Shown in Baby Yellow with Rosegold Details)



LIL OM (SHOWN IN SEAFOAM WITH LIME DETAILS)



LIL ROYALE (SHOWN IN WHITE WITH BLACK DETAILS)



LIL SOLO (SHOWN IN SKIN #2)

ALLELES[™] LiL Covers



LIL STEALTH (SHOWN IN LIME WITH PURPLE DETAILS)



LIL VICTORIAN (SHOWN IN SKIN #1 WITH ROSEGOLD DETAILS)



LIL WINGS (SHOWN IN RED WITH SILVER DETAILS)