


PROSTHETIC LINER SIZING CHARTS

HOW TO SIZE

Measurements taken
6cm from the distal end



AALPS

Tel +420 377 223 127 | Fax +420 377 223 010

www.easyliner.eu | info@easyliner.eu

Božkovské náměstí 17/21, 32600 Plzeň - Czech Republic

PART #	DESCRIPTION	THICKNESS (MM)	SIZE									
AKDT AKDTHD AKFR AKFRHD	Extreme Liner	3 or 6	16	20	24	26	28	32	38	44		
APDT-TP APFR-TP	Anterior Posterior Tapered Liner	6/3		20	24	26	28	32				
ECDTHD ECFRHD	ECO Liner	3 or 6		20	24	26	28	32	38	44		
ELDT ELFR	EasyLiner	3, 6, 6/3 or 9*	16	20	24*	26	28*	32	38	44		
GPDT GPFR	General Purpose	3 or 6	16	20	24	26	28	32	38	44		
SPDT SPDTHD SPFR SPFRHD	Superior Performance	3 or 6	16	20	24	26	28	32	38	44		
SADT SAFR SBFR	Smart Seal SADT	3 or 6		20	24	26	28	32	38			
	Smart Seal SAFR	3 or 6				26	28	32	38			
	Smart Seal SBFR	3 or 6		20	24	26	28					
EZF	EasyFlex	3 or 6	16	20	24	26	28	32	38	44		
VSDT	Liberty Liner	6		20	24	26	28	32	38			
PSI	Proximal Seal Interface	6 or 6/3	16	20	24	26	28	32	38	44		
CIRCUMFERENCE (CM)			16-19	20-23	24-25	26-27	28-31	32-37	38-44	44-53		

*9mm only available in sizes 24 and 28

PART #	DESCRIPTION	THICKNESS (MM)	SIZE										
LSC LSL	Silicone Pro	3	20	22	24	26	28	30	32	35	38	41	44
CIRCUMFERENCE (CM)			20-22	22-24	24-26	26-28	28-30	30-32	32-35	35-38	38-41	41-44	44-47

PART #	DESCRIPTION	THICKNESS (MM)	SIZE									
ELPX	EasyLiner SuperStretch	6	10	16	20	24	26	28	32	38		
CIRCUMFERENCE (CM)			10-15	16-19	20-23	24-25	26-27	28-31	32-37	38-43		