

COMPRESSION SOLUTIONS



CONTENTS

Section 1 - Juzo Compression Range

- Product characteristics 3
- Stockings 4
- Arm Sleeves 14
- Gloves & Toe Pieces 18
- Decongestive Therapy 22
- Compression Wraps 25
- Juzo Accessories 32



Section 2 - Jobst Compression Range

- Arm Sleeves & Gloves 36
- Lower Limb garments 43
- Jobst Farrowwraps 50
- Bellisse Bras 61
- Head & Neck 64
- Bandages & Accessories 65



Section 3 - Knitrite Compression Range

- Corespun Compression Socks 67
- Therafirm Compression Stockings 69
- Diabetic Crew Socks 71



Section 4 - LymphaPress Compression Pumps

- LymphaPress Advantages 74
- 8 Cell Pumps 75
- 12 Cell Pumps 77
- Garments 79



Section 5 - LymphaTouch Negative Pressure Device

- One Device Multiple Uses 92





VERSATILE FABRICS

For effective compression therapy



Medical compression stockings from Juzo guarantee correct pressure values and a well-defined pressure gradient. The pressure decreases continuously from the ankle section upwards towards the heart, and also towards the toes. We can provide a wide variety of garments, from below-knee stockings to pantyhose, from glove to face mask, depending on the therapeutic application.

Compression fabrics are categorized into 4 compression classes (CCL) based on the level of pressure. Possible indications can be, for example:

18 - 21 mmHg (CCL 1) // mild compression

- Minor obstructive disorders accompanied by feelings of tiredness or heaviness in the legs
- Isolated varicose vein formation
- Preventative measure during pregnancy
- Therapeutic treatment of scars

23 - 32 mmHg (CCL 2) // moderate compression

- Varicose veins
- Leg swelling during pregnancy
- After varicose vein sclerosis or operations
- For the long-term treatment of healed leg ulcers
- Propensity for oedema formation exists

34 - 46 mmHg (CCL 3) // strong compression

- Varicose veins accompanied by damage to the deep leg veins
- After healing of venous crural ulcers
- When a severe swelling tendency exists
- Chronic venous insufficiency
- Severe propensity for oedema formation

> 49 mmHg (CCL 4) // very strong compression

- Severe lymphatic obstruction
- Elephantiasis
- Lymphoedemas with fibrosis

All Juzo ranges can be custom made.

PRODUCT CHARACTERISTICS

Silicone borders



Top border

Top border without silicone - like the usual top of a sock. For fastening of thigh stockings (without silicone border) adhesive lotion or another attachment possibility may be necessary.



Elastic border

Durable border without silicone beads (dots). It stabilises the top border but has no adhesive properties and used for AD stockings and sleeves.



Silicone border

The silicone beads (dots) provide additional skin-friendly adhesion and prevent the stocking from slipping. The small beads allow the air to circulate well. The discreet silicone border is also suitable for male compression stocking wearers. For AD, CG and SG: thin silicone border (3,5 cm).



Decorative silicone border

The decorative silicone border has the same characteristics as the silicone border. This border option is especially demanded by women because there is a pretty but inconspicuous pattern on the outside of the border.



Lace border

A border option with a 3-row, continuous circular silicone coating on the inside. The silicone has excellent skin-friendly adhesion characteristics. Ladies appreciate the fine-stocking look and the elegant appearance due to the woven in floral* pattern.



Balance silicone border

The new balance silicone border with its micro-nub coating made from skin-friendly silicone has been specially tailored to the requirements of sensitive skin.

FINE FASHIONABLE STOCKINGS

JUZO INSPIRATION

Glossy with a fine-stocking look

The velvety transparency and exceptionally fine knit make this compression stocking a real trendsetter. Thanks to its fine-stocking look, Juzo Inspiration is perfectly suited for everyday use and is a versatile addition to your wardrobe. The fashion highlight of this range is the "Vision" pattern, which can be incorporated on request. The compression pantyhose with a special pregnancy body part optimally supports both the expectant mother and the baby.

- Fine-meshed fabric
- Available in Ccl.1, Ccl.2
- Article Number: 2701/2702



COLOURS

Juzo Inspiration Juzo Soft

Sugar	
Sesame	
Almond	
Cinnamon	
Cacao	
Poppy seed	
Blueberry	
Black pepper	

JUZO SOFT




Super soft for feel good moments

The fine microfibre fabric with its optimum moisture transport guarantees a comfortable skin environment and a pleasant, soft feeling, on the skin. The Juzo Soft combines refined fashion with ultimate in wearing comfort. You can make a fashion statement with the Juzo Soft leggings which can be combined with every outfit to ensure you feel your best!

- Soft fabric
- Available in Ccl.1, Ccl.2
- Article Number: 2001/2002



Juzo Inspiration

Trend Colours	
Batik Collection	
Dip Dye Collection	

Juzo Soft

Trend Colours	
Batik Collection	



Circular
Knitting

FINE FASHIONABLE STOCKINGS

SIZE CHARTS

Juzo Inspiration

Art no: 2701, 2702 Length in cm		Size: Circumferences in cm					
		I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
AG / AT / BT Normal 72.0 - 83.0 Short 66.0 - 72.0 Extra Short 58.0 - 66.0		Silicone border - Decorative border - Balance silicone border					
		Wide: Silicone border - Decorative silicone border - Balance silicone border					
		Pantyhose - Lace border - Top border					
AD Normal 40.0 - 46.0 Short 34.0 - 39.0 Extra Short 31.0 - 34.0		Thin silicone border - Top Border - Thin balance silicone border					
		Thin silicone border - Top Border - Thin balance silicone border					
		Thin silicone border - Top Border - Thin balance silicone border					
23.5 - 26.5 15.0 - 18.0		Thin silicone border - Top Border - Thin balance silicone border					
		Thin silicone border - Top Border - Thin balance silicone border					
Shoe Sizes: 37 - 42 (UK: 4 - 8)							

Juzo Soft

Art no: 2001, 2002 Length in cm		Size: Circumferences in cm					
		I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
AG / AT Normal 72.0 - 83.0 Short 66.0 - 72.0 Extra Short 58.0 - 66.0		Silicone border - Decorative silicone border - Balance silicone border					
		Wide: Silicone border - Decorative silicone border - Balance silicone border					
		Pantyhose - Leggings* - Lace border - Top border					
AD Normal 40.0 - 46.0 Short 34.0 - 39.0 Extra Short 31.0 - 34.0		Thin silicone border - Top Border - Thin balance silicone border					
		Thin silicone border - Top Border - Thin balance silicone border					
		Thin silicone border - Top Border - Thin balance silicone border					
23.5 - 26.5 15.0 - 18.0		Thin silicone border - Top Border - Thin balance silicone border					
		Thin silicone border - Top Border - Thin balance silicone border					
Shoe Sizes: 37 - 42 (UK: 4 - 8)		*Juzo Soft Leggings not in extra short					

INTENSIVE THERAPY STOCKINGS

JUZO MOVE

Pleasant wearing comfort and a high working pressure

The high working pressure offers optimal support for an active day to day life. The high microfibre content ensures the stocking is comfortable to wear while the short stretch fabric reduces slippage.

Thanks to it's functional zones, Juzo Move compression stockings provide real therapeutic relief for any patient and makes wearing more enjoyable and tolerable.


- Short stretch fabric
- Open or closed toe
- Available in Ccl.1, Ccl.2
- Range of borders available
- Article Number: 3611/3612



COLOURS

Juzo Move

Almond 


Black pepper 



Circular Knitting

SIZE CHARTS

Juzo Move

Art no: 3611, 3612 Length in cm		Size: Circumferences in cm					
		I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
AG Normal 72.0 - 83.0 Short 66.0 - 72.0 Extra Short 58.0 - 66.0		Silicone border - Decorative silicone border					
		43.0 - 49.0	47.0 - 53.0	49.0 - 57.0	55.0 - 63.0	59.0 - 67.0	64.0 - 72.0
AD Normal 40.0 - 46.0 Short 34.0 - 39.0 Extra Short 31.0 - 34.0		Wide: Silicone border - Decorative silicone border					
		50.0 - 56.0	54.0 - 62.0	58.0 - 66.0	64.0 - 73.0	68.0 - 77.0	73.0 - 82.0
		Lace border - Top border					
		43.0 - 56.0	47.0 - 62.0	49.0 - 66.0	55.0 - 73.0	59.0 - 77.0	64.0 - 82.0
		43.0 - 53.0	47.0 - 59.0	49.0 - 63.0	55.0 - 70.0	59.0 - 74.0	64.0 - 79.0
		28.0 - 34.0	30.0 - 37.0	33.0 - 41.0	35.0 - 43.0	41.0 - 51.0	44.0 - 55.0
		25.0 - 30.0	27.0 - 33.0	30.0 - 36.0	33.0 - 39.0	35.0 - 41.0	37.0 - 43.0
		27.0 - 34.0	29.0 - 37.0	33.0 - 41.0	36.0 - 45.0	39.0 - 47.0	42.0 - 50.0
		24.0 - 26.0	25.0 - 29.0	27.0 - 33.0	31.0 - 36.0	34.0 - 39.0	36.0 - 41.0
		18.0 - 20.0	20.0 - 22.0	22.0 - 24.5	24.5 - 27.0	27.0 - 30.0	30.0 - 33.0
		26.5 - 30.5	29.5 - 32.5	31.5 - 35.0	34.0 - 37.5	39.5 - 41.0	39.5 - 45.0
		15.0 - 18.0	17.5 - 22.0	19.5 - 24.0	21.5 - 26.5	24.0 - 29.0	26.5 - 33.0

Shoe Sizes: 37 - 42 (UK: 4 - 8)

INTENSIVE THERAPY STOCKINGS

JUZO ADVENTURE

The allrounder for active people

This stocking is perfectly tailored for the requirements of active people. With its attractive look, stable and durable structure, and opaque fabric, this stocking has been designed for use in a wide range of outdoor activities and sports. Thanks to the high cotton content of Juzo Adventure, it is also exceptionally skin friendly, making it comfortable to wear in the office as well.

- Sport
- Ribbed structure
- With cotton
- Closed toe
- Available in Ccl.1, Ccl.2
- Article Number: 3521/3522



COLOURS

Juzo Adventure

Salt-Pepper	
Almond	
Cacao	
Poppy seed	
Blueberry	
Black pepper	




Circular Knitting

SIZE CHARTS

Juzo Adventure

Art no: 3521, 3522 Length in cm	Size: Circumferences in cm					
	I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
eZ AD Normal 40.0 - 46.0 Short 34.0 - 39.0 Extra Short 31.0 - 34.0	cD 25.0 - 30.0	27.0 - 36.0	30.0 - 36.0	33.0 - 39.0	35.0 - 41.0	37.0 - 43.0
	cC 27.0 - 34.0	29.0 - 37.0	33.0 - 41.0	36.0 - 45.0	39.0 - 47.0	42.0 - 50.0
	cB' 24.0 - 26.0	25.0 - 29.0	27.0 - 33.0	31.0 - 36.0	34.0 - 39.0	36.0 - 41.0
	cB 18.0 - 20.0	20.0 - 22.0	22.0 - 24.5	24.5 - 27.0	27.0 - 30.0	30.0 - 33.0
	cY 26.5 - 30.5	29.5 - 32.5	31.5 - 35.0	34.0 - 37.5	36.5 - 41.0	39.5 - 45.0
eZ 25.0 - 29.0 (short 23.0 - 25.0)	cA 17.5 - 22.0	19.5 - 24.0	21.5 - 26.5	24.0 - 29.0	26.5 - 33.0	29.5 - 37.0

Shoe Sizes: Short Foot: 36 - 39 (UK: 3.5 - 5.5), Std. Foot: 39 - 45 (UK: 5.5 - 10)



INTENSIVE THERAPY STOCKINGS

JUZO DYNAMIC

Durability you can count on

Juzo Dynamic, with its tight and unobtrusive weave, is visually discreet. It offers a perfect fit, and retains its form even during high levels of physical exertion. This makes it the ideal circular knitted stocking for demanding requirements.

- Durable
- Available in Ccl.1, Ccl.2, Ccl.3
- Article Number: 3511/3512/3513





COLOURS

Juzo Dynamic

Sesame	
Almond	
Poppy seed	
Blueberry	
Black pepper	

Juzo Dynamic Cotton

Almond	
Black pepper	



Circular Knitting

JUZO DYNAMIC COTTON

With cotton for skin friendliness

This stocking, with comfortable cotton on the inside, is especially skin-friendly and feels pleasant on the skin even with skin irritations and allergies. The cotton, a natural fibre, optimally supports the functional characteristics of Juzo Dynamic.

- With cotton
- Durable
- Available in Ccl.1, Ccl.2
- Article Number: 3511/3512



INTENSIVE THERAPY STOCKINGS

SIZE CHARTS

Juzo Dynamic / Juzo Dynamic Cotton

Art no: 3511, 3512, 3513 / 3511, 3512 Cotton		Size: Circumferences in cm					
Length in cm		I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
øG AG / AT / BT Normal 72.0 - 83.0 Short 66.0 - 72.0 Extra Short 58.0 - 66.0		Silicone border - Decorative silicone border					
		43.0 - 49.0	47.0 - 53.0	49.0 - 57.0	55.0 - 63.0	59.0 - 67.0	64.0 - 72.0
		Wide: Silicone border - Decorative silicone border					
		50.0 - 56.0	54.0 - 62.0	58.0 - 66.0	64.0 - 73.0	68.0 - 77.0	73.0 - 82.0
		Pantyhose - Lace border - Top border					
		43.0 - 56.0	47.0 - 62.0	49.0 - 66.0	55.0 - 73.0	59.0 - 77.0	64.0 - 82.0
øZ AD Normal 40.0 - 46.0 Short 34.0 - 39.0 Extra Short 31.0 - 34.0		43.0 - 53.0	47.0 - 59.0	49.0 - 63.0	55.0 - 70.0	59.0 - 74.0	64.0 - 79.0
		28.0 - 34.0	30.0 - 37.0	33.0 - 41.0	35.0 - 43.0	41.0 - 51.0	44.0 - 55.0
		25.0 - 30.0	27.0 - 33.0	30.0 - 36.0	33.0 - 39.0	35.0 - 41.0	37.0 - 43.0
		27.0 - 34.0	29.0 - 37.0	33.0 - 41.0	36.0 - 45.0	39.0 - 47.0	42.0 - 50.0
		24.0 - 26.0	25.0 - 29.0	27.0 - 33.0	31.0 - 36.0	34.0 - 39.0	36.0 - 41.0
		18.0 - 20.0	20.0 - 22.0	22.0 - 24.5	24.5 - 27.0	27.0 - 30.0	30.0 - 33.0
øZ 23.5 - 26.5		26.5 - 30.5	29.5 - 32.5	31.5 - 35.0	34.0 - 37.5	39.5 - 41.0	39.5 - 45.0
øA 15.0 - 18.0		17.5 - 22.0	19.5 - 24.0	21.5 - 26.5	24.0 - 29.0	26.5 - 33.0	29.5 - 37.0

Shoe Sizes: 37 - 42 (UK: 4 - 8)

Have a question?

With an experienced team of allied health professionals right across Australia, we can assist with all of your clinical enquiries.

Have a question? Call us on [1300 672 937](tel:1300672937)



JUZO ULCER PRO

JUZO ULCER PRO

For the treatment of crural ulcers

Using its long-standing experience as a manufacturer of compression garments and working in close collaboration with medical technologists and therapists, Juzo has developed a compression stocking system that makes it easier to treat patients with ulcus cruris venosum. The treatment set comprises two types of compression stockings that are used in combination - one Ulcer Expert and two Ulcer Liners which hold wound dressings and pads securely in place. The compression stocking system has a lot of compelling features, amongst others it fits securely even during exercise and it can be worn with normal shoes.

- Stocking system
- Equivalent Ccl.3
- Article Number: 3073



COLOURS

Juzo Ulcer Liner

White



Juzo Ulcer Expert

Almond



Circular
Knitting
(Liner)



Flat
Knitting
(Expert)

JUZO ULCER LINER

The circular knitted Juzo Ulcer Liner stocking is suitable for long-term compression - 24 hours a day.

- Ccl.1
- 3 Liner Per Pack
- Article Number: 2571



JUZO ULCER EXPERT

The flat knitted Juzo Ulcer Expert cover stocking, enhances the compression effect during the active phase. It is worn over the liner stocking during the day to increase the level of compression.

- Ccl.2
- Article Number: 3022



JUZO ULCER PRO

SIZE CHARTS

Juzo Ulcer Pro

Art no: 3073 to 2571		Size: Circumferences in cm					
Length in cm		1 / XS	2 / S	3 / M	4 / L	5 / XL	6 / XXL
øZ AD		27.0 - 29.5	29.5 - 32.5	32.5 - 35.5	35.5 - 39.0	38.5 - 42.0	41.5 - 45.5
Normal	40.0 - 46.0						
Short	34.0 - 39.0						
		cD					
		cC					
		cB'					
		cB					
		cY					
øZ Juzo Ulcer Liner	23.5 - 26.5						
øZ Juzo Ulcer Expert	15.0 - 18.0						
		cA					

Shoe Sizes: 37 - 42 (UK: 4 - 8)

Stay up to Date

Like to stay in the know? Subscribe to the latest promotions, education and news from OPC Health!

To subscribe, visit opchealth.com.au



THE EXPERTS WITH SEAM

JUZO EXPERT

Pleasantly soft for tailored treatment

Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement and is very comfortable to wear. This special fabric is ideally suited for full-body garments and, thanks to its certified UV protection factor of 80 (18 - 32 mmHg, CCL 1 - 2), is also ideal for use in scar therapy.

- UV protection
- Flat seam
- Ccl.2 available in Standard Sizes with Top Border, open toe
- Ccl.2 available in Standard Sizes with Silicone Border, open toe
- MAX WIDTH available in Ccl.2 (borders limitations as above)
- Available in AD and AG
- Article Number: 3022



COLOURS

Juzo Expert Juzo Expert Cotton

Almond



Flat
Knitting

JUZO EXPERT COTTON

Natural fibres for extra softness

This stocking is made from almost 50 percent of cotton and covers the skin with ultrapure cotton. The natural fibres make it extremely soft and smooth. Juzo Expert Cotton is therefore especially comfortable to wear for sensitive, delicate skin.

- With cotton
- Flat seam
- Ccl.2 available in custom made.
- Article Number: 3022



THE EXPERTS WITH SEAM

SIZE CHARTS

Juzo Experts

Art no: 3022 Length in cm		Size: Circumferences in cm					
		I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
		Top border - Thin silicone border					
		29.5 - 32.5	32.5 - 35.5	35.5 - 39.0	38.5 - 42.0	41.5 - 45.5	45.5 - 48.5
ℓG	AD/AD max	cD					
Normal	40.0 - 46.0	27.0 - 29.5	29.5 - 32.5	32.5 - 35.5	35.5 - 39.0	38.5 - 42.0	41.5 - 45.5
Short	34.0 - 39.0	31.5 - 34.5	34.5 - 37.5	37.5 - 41.0	40.5 - 44.5	43.0 - 47.0	45.5 - 49.5
		cC					
		29.0 - 31.5	31.5 - 34.5	34.5 - 37.5	37.5 - 41.0	40.5 - 44.5	43.0 - 47.0
		cB'					
		26.0 - 28.5	28.5 - 31.0	30.5 - 33.5	33.0 - 36.5	35.5 - 38.0	37.5 - 41.0
		cB					
		23.5 - 25.5	25.5 - 28.0	27.5 - 30.0	30.0 - 33.0	32.5 - 35.5	35.0 - 38.5
ℓA		cY					
Normal	15.0 - 17.0	18.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0	26.0 - 28.5	28.0 - 30.5
Long	17.0 - 19.0	28.0 - 30.5	29.0 - 31.5	31.0 - 34.0	33.0 - 36.0	35.0 - 38.5	37.0 - 40.5
		cA					
		18.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0	26.0 - 28.5	28.0 - 30.5

Except max. sizes

Juzo Experts

Art no: 3022 Length in cm		Size: Circumferences in cm					
		I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
		Silicone border					
ℓG	AG	cG					
Normal	72.0 - 83.0	44.0 - 48.0	46.0 - 50.5	50.0 - 55.0	54.0 - 59.0	58.0 - 63.5	62.0 - 68.0
Short	66.0 - 72.0	cF					
		41.0 - 45.0	42.5 - 46.5	45.5 - 50.0	48.5 - 53.0	51.5 - 56.5	54.5 - 59.5
		cE					
		30.0 - 33.0	32.5 - 35.5	35.5 - 39.0	38.5 - 42.0	41.5 - 45.5	44.5 - 48.5
		cD					
		27.0 - 29.5	29.5 - 32.5	32.5 - 35.5	35.5 - 39.0	38.5 - 42.0	41.5 - 45.5
		cC					
		29.0 - 31.5	31.5 - 34.5	34.5 - 37.5	37.5 - 41.0	40.5 - 44.5	43.0 - 47.0
		cB'					
		23.5 - 25.5	25.5 - 28.0	27.5 - 30.0	30.0 - 33.0	32.5 - 35.5	35.0 - 38.5
		cB					
		18.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0	26.0 - 28.5	28.0 - 30.5
		cY					
		28.0 - 30.5	29.0 - 31.5	31.0 - 34.0	33.0 - 36.0	35.0 - 38.5	37.0 - 40.5
ℓA		cA					
AG Normal	17 - 19 cm	18.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0	26.0 - 28.5	28.0 - 30.5
AG Short	15 - 17 cm						



COMPRESSION SLEEVES

JUZO SOFT

Super soft for feel good moments











The fine microfibre fabric with its optimum moisture transport guarantees a comfortable skin environment and a pleasant, soft feeling on the skin. Juzo Soft combines refined fashion with the ultimate in wearing comfort so you can feel your best.

- Soft fabric
- Available either as a soft top border or with silicone
- Available in Ccl.1, Ccl.2
- Max width available
- Article Number: 2001/2002



COLOURS

Soft

Sugar	
Sesame	
Almond	
Cinnamon	
Cacao	
Poppy seed	
Blueberry	
Black pepper	
Trend Colours	
Batik Collection	

JUZO DYNAMIC

Durability you can count on - also with silver





Juzo Dynamic, with its tight and unobtrusive weave, is visually discreet. It offers a perfect fit, and retains its form even during high levels of physical exertion. This makes it the ideal circular knitted sleeve for demanding requirements.

The Juzo Dynamic Silver with a continuous knitted-in silver thread has a natural antibacterial effect. Price on application. Custom made.

- Durable
- With silver
- Available in Ccl.1, Ccl.2
- Article Number: 3511/3512



Dynamic

Sugar	
Sesame	
Almond	
Cinnamon	

JUZO EXPERT

Desired support with exceptional softness

Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement

- Ccl.2
- Please indicate Left or Right
- Seam on outside of arm, with overheight and thin silicone border
- Article Number: 3022



Circular Knitting

Juzo Expert

Almond	
--------	---




Flat Knitting


COMPRESSION SLEEVES

SIZE CHARTS


Juzo Soft Compression Sleeve

Art no: 2001, 2002 - CG, CH Length in cm	Size: Circumferences in cm						
	I / XS	II / S	III / M	IV / L	V / XL	VI / XXL	
 lCG 41.0 - 47.0 lCE 21.0 - 25.0	cG	24.0 - 26.0	26.0 - 28.5	28.0 - 30.5	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0
	max	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0	36.0 - 39.5	38.0 - 41.5	40.0 - 44.0
	cF	23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	max	28.0 - 30.5	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0	36.0 - 39.5	38.0 - 41.5
	cE	23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	max	25.0 - 27.5	27.0 - 29.5	29.0 - 31.5	30.5 - 33.5	32.0 - 35.0	33.5 - 36.5
	cD	23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	max	25.0 - 27.5	27.0 - 29.5	29.0 - 31.5	00.5 - 33.5	32.0 - 35.0	33.5 - 36.5
	cC	14.5 - 16.0	16.0 - 17.5	17.5 - 19.0	19.0 - 20.5	20.5 - 22.0	22.0 - 24.0

Juzo Dynamic Compression Sleeve

Art no: 3511, 3512 - CG, CH, Length in cm	Size: Circumferences in cm						
	I / XS	II / S	III / M	IV / L	V / XL	VI / XXL	
 lCG 41.0 - 47.0 lCE 21.0 - 25.0 lAC 7.5	cG	24.0 - 26.0	26.0 - 28.5	28.0 - 30.5	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0
	max	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0	36.0 - 39.5	38.0 - 41.5	40.0 - 44.0
	cF	23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	max	28.0 - 30.5	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0	36.0 - 39.5	38.0 - 41.5
	cE	23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	max	25.0 - 27.5	27.0 - 29.5	29.0 - 31.5	30.5 - 33.5	32.0 - 35.0	33.5 - 36.5
	cD	23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	max	25.0 - 27.5	27.0 - 29.5	29.0 - 31.5	30.5 - 33.5	32.0 - 35.0	33.5 - 36.5
	cC	14.5 - 16.0	16.0 - 17.5	17.5 - 19.0	19.0 - 20.5	20.5 - 22.0	22.0 - 24.0

Juzo Expert Compression Sleeve

Art no: 3022 - CG Length in cm	Size: Circumferences in cm				
	2 / S	3 / M	4 / L	5 / XL	
 lCG 42.0 - 46.0 lCE 22.0 - 24.0	cG	27.5 - 31.0	31.0 - 35.0	35.0 - 39.0	39.0 - 43.0
	cF	24.5 - 27.5	27.5 - 30.5	30.5 - 34.5	34.5 - 38.5
	cE	23.0 - 25.5	25.5 - 28.0	28.0 - 31.0	31.0 - 34.0
	cD	22.0 - 24.0	24.0 - 27.0	27.0 - 30.0	30.0 - 33.0
	cC'	17.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0
	cC	16.0 - 18.0	18.0 - 20.0	20.0 - 22.0	22.0 - 24.0

COMPRESSION SLEEVES



JUZO SOFT

Super soft for feel good moments

The fine microfibre fabric with its optimum moisture transport guarantees a comfortable skin environment and a pleasant, soft feeling on the skin. Juzo Soft combines refined fashion with the ultimate in wearing comfort so you can feel your best.

- Soft fabric
- 24 Size range option
- Silicone border as standard
- 15-20, 20-30, & 30-40mmHg available
- Matching gauntlet & gloves available



COLOURS

Soft

Almond	
White C06	
Black C10	
Chocolate C53	
Cinnamon C57	

JUZO PRINT SERIES



Floral Purple Henna



Floral Grey



Mosaic Henna



Tiger Jungle



Koi Flowers Tattoo



Bird Henna



Butterfly Henna



Floral Tan



Boho Bright

Please scan QR Code to view all the signature prints.



JUZO TREND COLOURS


Please scan QR Code to view all the latest trend colours.




COMPRESSION SLEEVES

SIZE CHARTS

Juzo Soft Compression Sleeve

Length in cm	Measuring Point	Sizes					
		I	II	III	IV	V	VI
C-G/C-H Regular C-G Length <17" 43cm C-G/C-H Long C-G Length >17" 43cm	 cG max	26.5 - 33	28.5 - 35	30.5 - 37.5	32.5 - 39.5	34 - 42	36 - 44
	cG	21.5 - 26	23.5 - 28.5	25 - 31	27 - 33	29 - 35	30.5 - 37.5
	cE max	22.5 - 27.5	24.5 - 29.5	26 - 32	27.5 - 33.5	29 - 35	30 - 37
	cE max	20.5 - 25	22 - 27	23.5 - 28.5	25 - 30	26 - 32	27.5 - 33.5
	cC	14 - 15.5	15.5 - 17	17 - 18.5	17 - 18.5	20 - 21.5	21.5 - 23

Juzo Soft Gauntlet

Length in cm	Measuring Point	Sizes			
		I	II	III	IV
	cA	17.5 - 19.5	19.5 - 21	21 - 23	23 - 25.5
	cC	14.5 - 16.5	18 - 20	17 - 18.5	20 - 22

Clinical Equipment

We can supply you with a wide range of clinical equipment. From treatment tables to massage chairs, trolleys to stools and everything in between.

Visit opchealth.com.au to view our entire range




COMPRESSION GLOVES & TOE PIECES

JUZO CLASSIC SEAMLESS GLOVE



Classic seamless CGA

- Balance Silicone Border



Classic seamless AGA

- Balance Silicone Border
- Knitted on gauntlet (with open thumb)




Classic seamless AGO

- Balance Silicone Border
- Knitted on glove (with open fingers & thumb)

Art no: 2301, 2302 Size 1 - 6		Model
	CGA	AGA AGO

Juzo Classic Seamless Compression Sleeve

Art no: 2301, 2302 - CGA, AGO, AGA	Size: Circumferences in cm						
	I / XS	II / S	III / M	IV / L	V / XL	VI / XXL	
	cG	25.0 - 29.0	28.5 - 32.5	31.5 - 36.0	35.0 - 39.5	38.5 - 43.0	42.0 - 45.5
	max	29.0 - 33.0	32.5 - 36.5	35.5 - 40.0	39.0 - 43.5	42.5 - 47.0	46.0 - 49.5
cF		24.5 - 27.0	26.5 - 30.0	29.0 - 33.0	32.0 - 36.0	35.0 - 38.0	37.0 - 40.0
	max	27.5 - 30.0	29.5 - 33.0	32.0 - 36.0	35.0 - 39.0	38.0 - 41.0	40.0 - 43.0
cE		22.5 - 25.0	24.5 - 27.5	26.5 - 30.0	29.0 - 32.0	31.0 - 34.0	33.0 - 36.0
	max	24.5 - 27.0	26.5 - 29.5	28.5 - 32.0	31.0 - 34.0	33.0 - 36.0	35.0 - 38.0
cD		22.0 - 23.5	23.0 - 26.0	25.0 - 28.0	27.0 - 30.0	29.0 - 32.0	31.0 - 34.0
	max	24.0 - 25.5	25.0 - 28.0	27.0 - 30.0	29.0 - 32.0	31.0 - 34.0	33.0 - 36.0
cC		13.5 - 15.0	15.0 - 16.5	16.5 - 18.0	18.0 - 20.0	20.0 - 22.0	22.0 - 24.0
	cB	16.5 - 18.0	18.0 - 19.5	19.5 - 21.0	21.0 - 23.0	23.0 - 25.5	25.5 - 28.0
cA	16.5 - 18.0	18.0 - 19.5	19.5 - 21.0	21.0 - 23.0	23.0 - 25.5	25.5 - 28.0	

Quality mark valid only for 2302 CGA

COMPRESSION GLOVES & TOE PIECES

JUZO CLASSIC SEAMLESS GLOVE & GAUNTLET

Comfortable thanks to its seamless fabric




Juzo Classic Seamless Gloves are setting new standards in oedema therapy. The innovative seamless fabric is breathable, thereby assuring a comfortable skin environment. Being soft and smooth as well, it guarantees an especially high wearing comfort and easy donning and doffing. The 3D comfort zone in the sensitive area between the thumb and index finger prevents the fabric from creasing and protects the skin from irritation. It is designed not to restrict the thumb movement in any way.

- Available Ccl.1 & Ccl.2
- Universal
- Article Number: 2301/2302



COLOURS

Classic Seamless

Sesame	
Almond	
Black pepper	

JUZO CLASSIC SEAMLESS TOE CAP

Seamless compression for lymphology

The lengthwise knitted structure of the toes ensures additional stability and easy donning of the toe part. The 3D comfort zone in the sensitive area between the big toe and the ball conforms perfectly to anatomical features, allowing for optimal fit and preventing the fabric from creasing. This protects the skin from irritation and guarantees a high wearing comfort.

- Available Ccl.1 & Ccl.2
- Please indicate either left or right
- Article Number: 2301/2302



JUZO EXPERT (seamless)




Pleasantly soft for tailored treatment

The anatomy of the hand is very complex, so medical compression garments for the hand need to be very sophisticated. Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement and is very comfortable to wear.

- Available Ccl.1 & Ccl.2
- Gauntlet available
- Glove with either closed or open fingers
- Article Number: 3021NL/3022NL



Expert Seamless

Sesame	
Almond	
Seasonal black	



Seamless Knitting




Flat Knitting


COMPRESSION GLOVES & TOE PIECES

SIZE CHARTS


Juzo Classic Seamless Gauntlet & Glove

Art no: 2301, 2302, HDA, HFA	Size - Measurements in cm					
	1 / XS	2 / S	3 / M	4 / L	5 / XL	6 / XXL
 cA	16.5 - 18.0	18.0 - 19.5	19.5 - 21.0	21.0 - 23.0	23.0 - 25.5	25.5 - 28.0
cC	13.5 - 15.0	15.0 - 16.5	16.5 - 18.0	18.0 - 20.0	20.0 - 22.0	22.0 - 24.0

Juzo Classic Seamless Toe Cap


Art no: 2301, 2302, ZEA	Size - Measurements in cm					
	1 / XS	2 / S	3 / M	4 / L	5 / XL	6 / XXL
 cA	18.0 - 20.0	20.0 - 22.5	22.5 - 25.0	25.0 - 27.5	27.5 - 30.0	30.0 - 32.5
cA	18.0 - 20.0	20.0 - 22.5	22.5 - 25.0	25.0 - 27.5	27.5 - 30.0	30.0 - 32.5

Juzo Expert Seamless & Expert Compression Gauntlet

Art no: 3021, 3022 - HDA Length in cm	Size: Circumferences in cm						
	1 / XS	2 / S	3 / M	4 / L	5 / XL	6 / XXL	
 l* ① 2.0	cX	6.1	6.5	6.9	7.4	7.8	8.2
cA		17.5	19.0	20.5	22.0	23.5	25.0
cB		14.5	16.0	17.5	19.0	20.5	22.0

*Length thumb stub

Juzo Expert Seamless & Expert Compression Gloves

Art no: 3021, 3022 NL - HFA Length in cm	Size: Circumferences in cm						
	1 / XS	2 / S	3 / M	4 / L	5 / XL	6 / XXL	
 l* ① 3.0 ② 4.0 ③ 5.0 ④ 4.0 ⑤ 3.0	cZ	6.1	6.5	6.9	7.4	7.8	8.2
cX		6.2	6.7	7.1	7.6	8.0	8.5
cA		6.0	6.3	6.7	7.1	7.5	7.9
cB		5.7	6.0	6.3	6.7	7.1	7.4
cC		5.3	5.5	5.9	6.3	6.7	7.0
cC		17.5	19.0	20.5	22.0	23.5	25.0

*Length finger stub

JUZO CLASSIC SEAMLESS

CT ALMOND



FT MAN ALMOND



FT WOMAN ALMOND












Smooth for incipient oedema

Juzo Classic Seamless products are the ideal solution for starting oedema therapy. They are especially suitable for treating incipient lymphoedema where the anatomical situation is still normal. The innovative knitting technology enables the seamless manufacture of compression garments.

SIZE CHARTS

Juzo Classic Seamless

Art no: 2302 - BT, CT, FT, FTM Length in cm	Size: Circumferences in cm					
	I / XS	II / S	III / M	IV / L	V / XL	VI / XXL
 cT	65.0 - 75.0	75.0 - 85.0	85.0 - 95.0	95.0 - 105.0	105.0 - 115.0	115.0 - 125.0
 cH	85.0 - 95.0	95.0 - 105.0	105.0 - 115.0	115.0 - 125.0	125.0 - 135.0	135.0 - 145.0
 cG	43.0 - 56.0	47.0 - 62.0	51.0 - 66.0	55.0 - 73.0	59.0 - 77.0	64.0 - 82.0
 cF	43.0 - 53.0	47.0 - 59.0	51.0 - 63.0	55.0 - 70.0	59.0 - 74.0	64.0 - 79.0
 cE	30.0 - 37.0	33.0 - 41.0	36.0 - 44.0	38.0 - 48.0	40.0 - 50.0	44.0 - 54.0
 cD	26.0 - 33.0	28.0 - 36.0	32.0 - 39.0	35.0 - 43.0	37.0 - 47.0	39.0 - 50.0
 cC	27.0 - 37.0	30.0 - 40.0	33.0 - 43.0	36.0 - 47.0	39.0 - 51.0	42.0 - 54.0
 cB'	22.5 - 28.5	25.0 - 31.0	27.5 - 33.5	30.0 - 37.0	33.0 - 40.5	36.0 - 43.5
 cB	18.5 - 20.0	19.0 - 22.0	21.0 - 24.5	23.5 - 27.0	26.0 - 30.0	29.0 - 33.0

JUZO PRODUCTS FOR DECONGESTIVE THERAPY

JUZO SOFTCOMPRESS

Special products for oedema therapy

Juzo SoftCompress Products are employed as an effective treatment aid for compression bandaging in oedema treatment. Using the comfortable bandages and accessories can significantly improve the treatment outcome. The vertically-stitched foam channels create differences in pressure, thereby promoting lymph drainage, while at the same time preventing soft tissue fibrosis. The Juzo SoftCompress products are especially skin-friendly and breathable. The material is washable up to 60 °C.

- Washable up to 60 °C
- Very breathable



JUZO SOFTCOMPRESS BANDAGES

Special products for decongestion

Bandaging aids for arm and leg which are easy to use and ideally suited to self-therapy. It can be worn under short-stretch bandages or in combination with Juzo ACS Light or Juzo Compression Wrap.

- Washable up to 60°C
- Very breathable



JUZO SOFTCOMPRESS PADS

Targeted increase in pressure

Juzo SoftCompress Pads are employed as an effective treatment aid for oedema. The vertically-stitched foam channels create differences in pressure, thereby promoting lymph drainage, while at the same time preventing soft tissue fibrosis.

- Washable up to 60°C
- Very breathable



JUZO SOFTCOMPRESS COMPRESSION ROLL

Additional pressure

The thin, very soft Juzo SoftCompress Compression Rolls already produce compression pressure if applied suitably. The compression rolls distribute the pressure of the additionally applied short-stretch bandages, thereby avoiding constriction.

- Washable up to 60°C
- Very breathable



JUZO SOFTCOMPRESS

SIZE CHARTS

Juzo SoftCompress Bandage Arm

Art no: J6822 Length in cm	Size: Circumferences in cm					
	1 / XS	2 / S	3 / M	4 / L	5 / XL	6 / XXL
ℓG 46.0 - 51.0	cG 30.0 - 33.0	32.0 - 35.0	34.0 - 37.0	36.0 - 39.5	38.0 - 41.5	40.0 - 44.0
	cF 28.0 - 30.5	30.0 - 33.0	32.0 - 35.0	34.0 - 37.0	36.0 - 39.5	38.0 - 41.5
ℓE 21.0 - 25.0	cE 25.0 - 27.5	27.0 - 29.5	29.0 - 31.5	30.5 - 33.5	32.0 - 35.5	33.5 - 36.5
	cD 23.0 - 25.0	24.5 - 27.0	26.0 - 28.5	27.5 - 30.0	29.0 - 31.5	30.5 - 33.5
	cC 14.5 - 16.0	16.0 - 17.5	17.5 - 19.0	19.0 - 20.5	20.5 - 22.0	22.0 - 24.0

Juzo SoftCompress Bandage Lower Leg

Art no: J6830 Length in cm	Size: Circumferences in cm				
	1 / XS	2 / S	3 / M	4 / L	5 / XL
ℓD 38.0 - 46.0	cD 27.0 - 29.5	29.5 - 32.5	32.5 - 35.5	35.5 - 39.0	38.5 - 42.0
	cC 29.0 - 31.5	31.5 - 34.5	34.5 - 37.5	37.5 - 41.0	40.5 - 44.5
	cB' 23.5 - 25.5	25.5 - 28.0	27.5 - 30.0	30.0 - 33.0	32.5 - 35.5
	cB 18.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0	26.0 - 28.5
	cY 28.0 - 30.5	29.0 - 31.5	31.0 - 34.5.0	33.0 - 36.0	35.0 - 38.5
	cA 18.5 - 20.0	20.0 - 22.0	22.0 - 24.0	24.0 - 26.0	26.0 - 28.5

Practice Supplies

We can supply you with a wide range of clinical supplies. From wipes to gloves, sports tape to massage creams and everything in between!


Visit opchealth.com.au to view our entire range




SPECIAL PRODUCTS FOR OEDEMA THERAPY

SIZE CHARTS

Juzo SoftCompress Breast Pad

Art no: J6837	Size	Cup
	1	A - B
	2	C - D
	3	E - F


Juzo SoftCompress Compression Sheet

Art no: J6836	Size	Measurements in cm (breadth x height)
	1	65.0 x 65.0
	2	125.0 x 100.0


Juzo SoftCompress Breast bandage

Art no: J6838	Size	Cup	Breast Circumference in cm	Underbreast circumference in cm
	1	A	90 - 106.0	84.0 - 100.0
	2	B	90 - 106.0	84.0 - 100.0
	3	C	95 - 112.0	85.0 - 102.0
	4	D	95 - 112.0	85.0 - 102.0
	5	E	100 - 116.0	88.0 - 104.0
	6	F	100 - 116.0	88.0 - 104.0


Juzo SoftCompress Genital Pad Women

Art no: J6839	Size	Measurements in cm (breadth x height)
	1	20.0 x 20.0
	2	30.0 x 28.0
	3	40.0 x 32.0

Juzo SoftCompress Genital Pad Men

Art no: J6840	Size	Measurements in cm (breadth x height)
	1	20.0 x 28.0
	2	23.0 x 30.0
	3	26.0 x 32.0

Juzo SoftCompress Compression Roll

Art no: J6835	Size	Measurements in cm (breadth x height)
	1	400.0 x 8.0
	2	400.0 x 10.0
	3	400.0 x 12.0
	4	400.0 x 16.0

Thickness in mm: 5mm

JUZO COMPRESSION WRAP

JUZO COMPRESSION WRAP

For the treatment of *ulcus curis venosum*

The Juzo Compression Wrap is a short-stretch garment that can be employed as an alternative to conventional wrap-around bandages. Thanks to its short-stretch properties, it delivers a very high working pressure and a low resting pressure. The velcro fastener makes it easy to apply, and easy to adjust for fluctuations in body circumference.

The special indentations prevent "window oedema" by ensuring that no gaps can form between the straps. The Juzo Compression Wrap is reversible, with black on one side and beige on the other, to suit the individual patient's preference.

- Very high working pressure
- Reversible
- Arm and leg wraps available
- Working Pressure 60mmHg
- Washable at 40°C



COLOURS

Soft


Almond / Black



JUZO COMPRESSION WRAP


SIZE CHARTS

Juzo Compression Wrap - Foot


Art no: 6000 LF/LFL	Size - Measurements in cm			
	2 / S	3 / M	4 / L	5 / XL
	15.0 - 23.0	20.0 - 28.0	25.0 - 33.0	30.0 - 40.0
	19.0 - 24.0	24.0 - 29.0	29.0 - 34.0	34.0 - 40.0
normal	16	18	19.5	19.5
long	19	21	22.5	23

* 6cm proximal to the little toe Metatarsophalangeal joint
 ** The circumference can be adjusted using the extended tab at "cB" if required.


Juzo Compression Calf Wrap

Art no: 6000		Size: Circumferences in cm					
		2 / S	3 / M	4 / L	5 / XL	6 / XXL	7 / XXXL
	AD						
Normal	< 40.0	29.0 - 35.0	34.0 - 40.0	37.0 - 44.0	41.0 - 48.0	46.0 - 53.0	50.0 - 60.0
Short	> 40.0						
	max	35.0 - 41.0	40.0 - 46.0	44.0 - 51.0	48.0 - 55.0	53.0 - 60.0	60.0 - 70.0
		18.0 - 21.0	21.0 - 24.0	24.0 - 27.0	27.0 - 31.0	31.0 - 35.0	36.0 - 47.0

Juzo Compression Knee Wrap

Art no: 6000 LK	Circumference c	Standard Sizes - Circumferences in cm					
		2 / S	3 / M	4 / L	5 / XL	6 / XXL	7 / XXXL
		30.0 - 36.0	36.0 - 42.0	42.0 - 48.0	48.0 - 55.0	55.0 - 62.0	62.0 - 70.0

Juzo Compression Thigh Wrap

Art no: 6000 LT	Circumference c Length l	Size - Measurements in cm				
		2 / S	3 / M	4 / L	5 / XL	6 / XXL
	largest point of thigh	43.0 - 70.0	54.0 - 81.0	65.0 - 92.0	75.0 - 102.0	85.0 - 112.0
	10 cm above patella	35.0 - 55.0	48.0 - 70.0	60.0 - 82.0	70.0 - 92.0	80.0 - 102.0

JUZO COMPRESSION WRAP

Juzo Compression Wrap - Hand

Art no: 6000 UHR / UHL	Size - Circumferences in cm				
	1 / XS	2 / S	3 / M	4 / L	5 / XL
	15.5 - 17.5	17.5 - 19.5	19.5 - 22.0	22.0 - 24.5	24.5 - 28.0
	13.0 - 17.0	15.0 - 19.0	17.0 - 21.5	20.0 - 24.0	21.0 - 25.0

Juzo Compression Arm Wrap

Art. no. 6000 UH	Circumference c Length ℓ	Size - Measurements in cm				
		1 / XS	2 / S	3 / M	4 / L	5 / XL
	Upper arm	27.0 - 36.0	30.0 - 39.0	33.0 - 44.0	35.0 - 49.0	40.0 - 56.0
	Elbow	23.0 - 32.0	24.0 - 34.0	26.0 - 36.0	29.0 - 42.0	31.0 - 46.0
	Wrist	16.0 - 19.0	17.0 - 20.0	18.0 - 22.0	20.0 - 26.0	21.0 - 29.0

Have a question?

With an experienced team of allied health professionals right across Australia, we can assist with all of your clinical enquiries.

Have a question? Call us on [1300 672 937](tel:1300672937)



JUZO ADJUSTABLE COMPRESSION SYSTEM

JUZO ACS LIGHT

Efficient and effective adjustable compression system

The short-stretch compression wraps can be employed as an efficient supplementary treatment during the decongestion phase. The Juzo ACS light products can be adjusted for the individual wearer by means of velcro fasteners. Thanks to the incorporated lining, the wrap is easy to use, can be applied in next to no time and thus increases patient compliance.

- Easy to don
- High working pressure
- Durable
- Working Pressure 40mmHg
- Washable at 40°C
- Easily adjustable
- Available in Regular & Long versions



COLOURS

ACS light

Black



Almond



JUZO ADJUSTABLE COMPRESSION SYSTEM

SIZE CHARTS

Juzo ACS light Foot Wrap

Art no: J701LI	Circumference Length ℓ	Size - Measurements in cm				
		1 / S	2 / M	3 / L	4 / L	5 / XL
	Over forefoot	22.0 - 24.0	25.0 - 27.0	28.0 - 30.0	31.0 - 34.0	35.0 - 40.0
	Normal	16.5	18.0	19.5	21.0	23.0
	Long	18.5	20.0	21.5	23.0	25.0

Juzo ACS light Calf Wrap

Art no: J702LI	Circumference Length ℓ	Size - Measurements in cm				
		1 / S	2 / M	3 / L	4 / L	5 / XL
	Largest point of calf	35.0 - 45.0	40.0 - 52.0	48.0 - 60.0	55.0 - 65.0	60.0 - 70.0
	Ankle	22.0 - 27.0	25.0 - 32.0	30.8 - 38.0	37.0 - 44.0	45.0 - 52.0
	Normal	30.0	30.0	35.0	35.0	37.0
	Long	35.0	35.0	40.0	40.0	42.0

* below patella

Juzo ACS light Knee Wrap

Art no: J703LI	Circumference c	Size - Measurements in cm				
		1 / S	2 / M	3 / L	4 / L	5 / XL
	10cm above patella	40.0 - 55.0	50.0 - 65.0	60.0 - 75.0	70.0 - 85.0	85.0 - 100.0
	Below patella	35.0 - 45.0	40.0 - 52.0	48.0 - 60.0	55.0 - 65.0	60.0 - 70.0


Juzo ACS light Thigh Wrap

Art no: J704LI	Circumference c Length ℓ	Size - Measurements in cm				
		1 / XS	2 / S	3 / M	4 / L	5 / XL
	largest point of thigh	50.0 - 65.0	60.0 - 65.0	70.0 - 80.0	75.0 - 90.0	85.0 - 100.0
	10 cm above patella	45.0 - 52.0	53.0 - 59.0	60.0 - 66.0	67.0 - 73.0	74.0 - 78.0.0
	Normal	23.0	23.0	23.0	23.0	23.0
	Long	29.0	29.0	29.0	29.0	29.0

JUZO ADJUSTABLE COMPRESSION SYSTEM

SIZE CHARTS

Juzo ACS light Hand Wrap

Art no: J705LI	Circumference c	Size - Measurements in cm			
		1 / XS	2 / S	3 / M	4 / L
 cA	Over Knuckles	< 18.0	18.0 - 20.0	20.0 - 22.0	> 22.0


Juzo ACS light Arm Wrap

Art. no. J706LI	Circumference c Length ℓ	Size - Measurements in cm				
		1 / XS	2 / S	3 / M	4 / L	5 / XL
 cG cE cC ℓCG	Upper arm	25.0 - 28.0	29.0 - 33.0	34.0 - 39.0	40.0 - 45.0	46.0 - 50.0
	Elbow	23.0 - 26.0	26.0 - 30.0	29.0 - 34.0	33.0 - 39.0	38.0 - 42.0
	Wrist	14.0 - 16.0	17.0 - 19.0	20.0 - 22.0	23.0 - 25.0	24.0 - 27.0
	Normal	42.0 - 47.0	42.0 - 47.0	42.0 - 22.0	47.0 - 25.0	42.0 - 47.0
	Long	46.0 - 50.0	46.0 - 50.0	46.0 - 50.0	46.0 - 50.0	46.0 - 50.0

Stay up to Date

Like to stay in the know? Subscribe to the latest promotions, education and news from OPC Health!

To subscribe, visit [opchealth.com.au](https://www.opchealth.com.au)



JUZO LINER STOCKINGS SILVER

JUZO SILVER LINER STOCKINGS

With antibacterial effect, versatile

These soft-fabric stockings are especially comfortable to wear and are also suitable as a compression Aid for patients with extremely sensitive skin. The stockings offer permanent skin protection during slipping-on and taking-off and, thanks to the knitted-in silver thread, also have an antibacterial action. Juzo Silver Liner Stockings are ideally suited for incorporating wound dressings or securing special pressure pads.



COLOURS

Soft

Light Grey



SIZE CHARTS

Juzo Silver Liner Stockings

Art. nr.	Sizes	Circumferences in cm
J9640/1	S/M	Calf circumferences > 50 cm
J9640/2	I/XL	Calf circumferences < 50 cm

ARION SLIP-ON AND TAKING-OFF AIDS

SPECIAL GLOVES

High-quality textile glove with rubberised palm.

Long-life gloves for easier donning of compression garments.

- Code: J9301
- S, M, L or XL



ARION EASY-SLIDE ARM - CHILDREN

Donning aid for compression sleeves and stockings / pantyhose with open toes.

Smooth, two-layer material, very easy to use.

- Code: J1712
- XS or S



ARION SLIM-SLIDE

Donning and doffing aid for compression stockings / pantyhose.

Smooth, two-layer material.

Simply fasten to the ankle using the elastic strap with open toes.

- Code: J1705
- S, M, L or XL



ARION EASY-SLIDE LEG

Donning aid for compression stockings / pantyhose with open toes.

Smooth, two-layer material, very easy to use.

- Code: J1702
- S, M, L or XL



ARION EASY-OFF

Donning and doffing aid for compression stockings / pantyhose.

Smooth, two-layer material.

With velcro fastener for easy fixation above the ankle with open toes.

- Code: J1711
- S, M or L



ARION MAGNIDE

Donning aid for compression stockings / pantyhose with closed toes.

Smooth, two-layer material, very easy to use.

With innovative magnetic closure.

- Code: J1703
- M, L or XL



ARION EASY-SLIDE ARM

Donning aid for compression sleeves and elbow supports.

Smooth, two-layer material, very easy to use.

- Code: J1701
- M or L



ARION MAGNIDE 2-IN-1

Donning and doffing aid for compression stockings / pantyhose with closed toes.

Smooth, two-layer material, very easy to use.

With innovative magnetic closure.

- Code: J1709
- M, L or XL



JUZO PRESSURE AND LYMPH PADS

PRESSURE PADS BY DR. ROTTER

Covered with textile or silver-containing (antibacterial) material.

For DIY fitting.

- Code: J9400



JUZO LYMPHPAD SQUARE

Treatment aid for lymph- and lipoedema.

Synthetic latex, polyamide fleece on the rear.

Size 12 x 16cm (can be cut to any size).

- Code: J4311



PRESSURE PAD FOR THE HAND

Covered with textile or silver-containing (antibacterial) material.

To be inserted between back of the hand and glove.

- Code: J9420



JUZO LYMPH LOTION

Skin care lotion with 4% urea for daily use.

- Code: J9860



JUZO LYMPHPAD LINE

Treatment aid for lymph- and lipoedema.

Synthetic latex, polyamide fleece on the rear.

Size 12 x 16cm (can be cut to any size).

- Code: J9430



JUZO SPECIAL DETERGENT

For medical compression garments, all types of elastic fabrics and textiles, and high quality fabrics made from wool and silk, highly concentrated and economical to use.

- Code: J9000-250



JUZO ADHESIVE LOTION

Special roll-on adhesive for comfortable wearing of compression stockings and sleeves.

Size 60ml

- Code: J9800





COMPRESSION

To uplift your life



JOBST are passionate about providing compression therapy solutions that empower you to confidently live the life you want, without compromise - despite your lymphatic or venous condition. Everyone has different needs, so Jobst design their products to support those differences. Jobst's knowledge, insights and years of experience mean that they keep providing innovative products and services that suit your lifestyle, giving you the comfort, support and peace of mind that you truly deserve while wearing compression.

Jobst provide a wide variety of garments, from below-knee stockings to pantyhose and sleeves to gloves.

Compression fabrics are categorised into 3 compression classes (CCL) based on the level of pressure. Possible indications can be, for example:

15-20 mmHg (CCL1) mild compression

- Tired aching legs
- Minor ankle, leg and foot swelling
- Minor hand/arm swelling
- Minor varicose veins during pregnancy
- Long distance travel

20-30 mmHg (CCL2) moderate compression

- Tired aching legs
- Moderate to severe varicose veins during pregnancy
- Post vein surgery/post sclerotherapy
- Moderate tendency to oedema
- Superficial thrombophlebitis
- Helps prevent recurrence of ulcerations
- Upper limb lymphatic oedema

30-40 mmHg (CCL3) strong compression

- Severe varicose veins
- Lymphatic oedema
- Post vein surgery/post sclerotherapy
- Venous ulcers
- Manages manifestations of post thrombotic syndrome (PTS), chronic venous insufficiency CVI/CVD
- Aids in prevention of thrombotic syndrome

Jobst also offer a large range of custom made garments in a variety of fabrics.

JOBST® Bella™ LITE

JOBST BELLA LITE ARMSLEEVE AND GAUNTLET



Early detection and effective management of lymphoedema can reduce the symptom severity and improve quality of life.


The JOBST Bella Lite range is a family of products which are soft, air permeable, light weight and easy to don.


Available in 20-30mmHg, they provide an effective means for delivering compression in the earlier stages of oedema management.


The JOBST Bella Lite Armsleeve and Gauntlet can be worn together or separately to help maintain the limb size of lymphoedema patients.

SIZE CHARTS


JOBST Bella Lite Armsleeve


Wrist		S	M	L
	Long	75276-00	75277-00	75278-00
	Regular	75276-05	75277-05	75278-05
		14cm - 20cm	16cm - 22cm	18cm - 24cm

Mid-lower		S	M	L
	Long	75276-00	75277-00	75278-00
	Regular	75276-05	75277-05	75278-05
		18cm - 23cm	23cm - 29cm	29cm - 34cm

Mid-upper		S	M	L
	Long	75276-00	75277-00	75278-00
	Regular	75276-05	75277-05	75278-05
		22cm - 40cm	27cm - 45cm	33cm - 51cm

JOBST Bella Lite Gauntlet

Palm		S	M	L
		75282-00	75283-00	75284-00
		15cm - 19cm	19cm - 22cm	21cm - 26cm

Wrist		S	M	L
		75282-00	75283-00	75284-00
		14cm - 20cm	16cm - 22cm	18cm - 24cm

JOBST® Bella™ LITE


JOBST BELLA LITE COMBINED ARMSLEEVE AND GAUNTLET



JOBST Bella Lite Combined Sleeve provides the comfort and effectiveness of the Bella Lite armsleeve and gauntlet in one garment.

SIZE CHARTS


Lengths



Regular (up to 43 cm)
Long (over 43 cm)

JOBST® Bella™ Lite combination sleeve REGULAR

Palm	S	M	L
	77680-00	77680-01	77680-02
	15cm - 19cm	19cm - 22cm	22cm - 26cm




JOBST® Bella™ Lite combination sleeve LONG

S	M	L
77680-03	77680-04	77680-05
15cm - 19cm	19cm - 22cm	22cm - 26cm

JOBST® Bella™ Lite combination sleeve REGULAR

Mid-lower	S	M	L
	77680-00	77680-01	77680-02
	18cm - 23cm	23cm - 29cm	29cm - 34cm




JOBST® Bella™ Lite combination sleeve LONG

S	M	L
77680-03	77680-04	77680-05
18cm - 23cm	23cm - 29cm	29cm - 34cm

JOBST® Bella™ Lite combination sleeve REGULAR

Mid-lower	S	M	L
	77680-00	77680-01	77680-02
	22cm - 40cm	27cm - 45cm	33cm - 51cm



JOBST® Bella™ Lite combination sleeve LONG

S	M	L
77680-03	77680-04	77680-05
22cm - 40cm	27cm - 45cm	33cm - 51cm

JOBST® Bella™ LITE

JOBST BELLA LITE GLOVE




The JOBST Bella Lite Glove goes hand in hand with the JOBST Bella Lite sleeve.

Durable and latex free, the glove provides gentle compression for patients with mild to moderate oedema. With an ultra-thin fabric and fine seams, the glove provides great wearing comfort. The eight sizes provides a broad range of available options.

It is suitable for either right or left hand for added convenience in size selection. Fingers can also be cut to length without fraying.

MEASUREMENTS TABLE



Circumference in cm

Size	1	2	3	4
A	16.6-18.0 cm	16.6-18.0 cm	18.0-19.5 cm	18.0-19.5 cm
C	15.6-16.7 cm	16.7-18.1 cm	15.9-17.0 cm	17.0-18.5 cm
Code	71026-00	71026-01	71026-02	71026-03

Circumference in cm

Size	5	6	7	8
A	18.0-19.5 cm	19.5-21.5 cm	19.5-21.5 cm	21.5-23.1 cm
C	18.5-20.4cm	17.2-18.7 cm	18.7-20.6 cm	18.3-20.1 cm
Code	71026-04	71026-05	71026-06	71026-07

JOBST® Bella™ STRONG

JOBST BELLA STRONG ARMSLEEVE



Using the measurements from over 36,000 lymphoedema patients, the JOBST Bella Strong range has been able to develop an innovative sizing system. As a result, the armsleeves can fit approximately 80% of patients, and gauntlets are able to fit approximately 90% of patients.


Lymphoedema can be one of the most troubling complications after breast cancer treatment. Having access to optimal fitting compression garments can provide true confidence and enhance quality of life.

Available in 20-30mmHg, JOBST Bella Strong utilises advanced comfort technology to regulate humidity, keeping you cool and comfortable. Suitable for mild to moderate lymphoedema.

The JOBST Bella Strong Armsleeve and Gauntlet can be worn together or separately to help maintain the limb size of lymphoedema patients.

A silicone band on the armsleeve ensure a secure hold, preventing slippage and curling.

JOBST BELLA™ STRONG ARMSLEEVE MEASURING AND ORDERING GUIDE



All sizes are available in 2 lengths:
Regular up to 43 cm
Long greater than 43 cm

Sizing instructions

Finding the right size for each patient is easy:

1. Measure mid lower circumference (D). This leaves 2 or 3 choices
2. Measure axilla circumference (G). This leaves 1 or 2 choices
3. Measure wrist circumference (C). This leaves 1 or 2 choices

If there are still 2 options left, we recommend choosing the larger size for ease of donning, overall comfort and ultimate compliance. Either size selection will deliver the desired therapy.

Size	Wrist (C)	Mid Forearm (D)	Axilla (G)
1	14-18 cm	17-21 cm	22-31 cm
2	14-18 cm	17-21 cm	31-41 cm
3	14-18 cm	21-25 cm	25-34 cm
4	16-19 cm	21-25 cm	28-38 cm
5	17-20 cm	21-25 cm	34-43 cm
6	14-18 cm	25-29 cm	29-38 cm
7	17-20 cm	25-29 cm	33-43 cm
8	18-21 cm	25-29 cm	38-48 cm
9	16-19 cm	29-33 cm	32-42 cm
10	19-22 cm	29-33 cm	36-46 cm

20-30 SILICONE NATURAL REGULAR

Size	Code	Size	Code
1	72736-01	6	72736-06
2	72736-02	7	72736-07
3	72736-03	8	72736-08
4	72736-04	9	72736-09
5	72736-05	10	72736-10

20-30 SILICONE NATURAL LONG

Size	Code	Size	Code
1	72737-01	6	72737-06
2	72737-02	7	72737-07
3	72737-03	8	72737-08
4	72737-04	9	72737-09
5	72737-05	10	72737-10

JOBST® Bella™ STRONG

JOBST BELLA STRONG GAUNTLET



Using the measurements from over 36,000 lymphoedema patients, the JOBST Bella Strong range has been able to develop an innovative sizing system. As a result, the armsleeves can fit approximately 80% of patients, and gauntlets are able to fit approximately 90% of patients.

Lymphoedema can be one of the most troubling complications after breast cancer treatment. Having access to optimal fitting compression garments can provide true confidence and enhance quality of life.

Available in 20-30mmHg, JOBST Bella Strong utilises advanced comfort technology to regulate humidity, keeping you cool and comfortable. Suitable for mild to moderate lymphoedema.

The JOBST Bella Strong Armsleeve and Gauntlet can be worn together or separately to help maintain the limb size of lymphoedema patients.

A silicone band on the armsleeve ensure a secure hold, preventing slippage and curling.

JOBST BELLA™ STRONG GAUNTLET MEASURING GUIDE - 2 EASY STEPS


20-30 NATURAL

Sizing instructions

1. Measure palm circumference, this leaves 1 or 2 choices
2. Measure wrist circumference (select gauntlet)

Size	Palm (B)	Wrist (C)	Code
1	15-17 cm	14-18 cm	73694-40
2	17-19 cm	14-18 cm	73694-41
3	17-19 cm	18-21 cm	73694-42
4	19-21 cm	14-18 cm	73694-43
5	19-21 cm	18-21 cm	73694-44
6	21-23 cm	16-19 cm	73694-45
7	21-23 cm	19-22 cm	73694-46
8	23-26 cm	19-22 cm	73694-47

JOBST BELLA STRONG GLOVE



Leveraging the advanced flat knitting technology only ever seen in our custom garments, patients can now access that near custom feel in a ready to wear product.

The glove is available in eight different sizes, covering a broad range of patients. It is a perfect complement to the JOBST Bella Strong Armsleeve.

MEASUREMENTS TABLE

Circumference in cm

Size	1	2	3	4
A	16.6-18.0 cm	16.6-18.0 cm	18.0-19.5 cm	18.0-19.5 cm
C	15.6-16.7 cm	16.7-18.1 cm	15.9-17.0 cm	17.0-18.5 cm
Code	75291-10	75291-11	75291-12	75291-13
Size	5	6	7	8
A	18.0-19.5 cm	19.5-21.5 cm	19.5-21.5 cm	21.5-23.1 cm
C	18.5-20.4 cm	17.2-18.7 cm	18.7-20.6 cm	18.3-20.1 cm
Code	75291-14	75291-15	75291-16	75291-17

JOBST® ELVAREX

ELVAREX® READY-TO-WEAR SLEEVES




Take advantage of the trusted flat knit compression with Jobst Elvarex.

Jobst Elvarex RTW gloves, gauntlets and sleeves, offer effective oedema management for those requiring firm support. Available in 20-30 mmHg.


For those patients who do not fit into the sizing parameters, please see our extensive custom-made offering.

Elvarex Ready-to-Wear sleeves are available in 3 sizes and have a compression level of 20-30 mmHg.

Styles available:



Bias finish (no strap)



Bra loop




Body strap

ARMSLEEVE LENGTH SUITABLE FOR UP TO 42CM FOR CG AND 45CM FOR CG WITH ATTACHMENT

Measurement Point	Small	Medium	Large
C (Wrist)	16 - 18 cm	18 - 20 cm	20 - 23 cm
D (Forearm)	23 - 26 cm	26 - 28 cm	28 - 30 cm
E (Elbow)	24 - 26 cm	26 - 29 cm	29 - 32 cm
F (Mid Upper arm)	26 - 29 cm	29 - 33 cm	33 - 37 cm
G (Axilla)	31 - 34 cm	34 - 38 cm	38 - 44 cm

Sleeve with bias		Sleeve with bra loop	
Small Right	76101-00	Small Right	79435-00
Medium Right	76101-02	Medium Right	79435-02
Large Right	76101-04	Large Right	79435-04
Small Left	76101-01	Small Left	79435-01
Medium Left	76101-03	Medium Left	79435-03
Large Left	76101-05	Large Left	79435-05

Sleeve with body strap			
Small Left	79435-07	Small Right	79435-06
Medium Left	79435-09	Medium Right	79435-08
Large Left	79435-11	Large Right	79435-10



JOBST® ELVAREX


ELVAREX® READY-TO-WEAR GLOVES & GAUNTLETS




As with sleeves, Ready-to-Wear gloves and gauntlets are available in three sizes: small, medium and large. The following diagram illustrates the measuring points.

SIZE CHARTS

Style available:



Gauntlets




Glove

LENGTHS A-B 3CM, A-C 9CM, A-C1 12CM

Measurement Point	Small	Medium	Large
A (Metacarpal heads)	16 - 18 cm	17.5 - 19.5 cm	19.5 - 20.5 cm
B (Widest Palm)	17 - 19 cm	19 - 21 cm	21 - 23 cm
C (Wrist)	16 - 18 cm	18 - 20 cm	20 - 22 cm
Thumb base	5.0 - 6.0 cm	6.0 - 7.0 cm	7.0 - 8.0 cm
Fingers 2 base	5.0 - 6.0 cm	6.0 - 7.0 cm	6.0 - 7.0 cm
Fingers 3 base	5.0 - 6.0 cm	6.0 - 7.0 cm	6.5 - 7.0 cm
Fingers 4 base	4.5 - 5.5 cm	5.5 - 6.5 cm	6.5 - 7.5 cm
Fingers 5 base	4.0 - 5.0 cm	5.0 - 6.0 cm	6.0 - 7.0 cm

Gauntlet		Glove	
Small Right	76091-00	Small Right	76091-06
Medium Right	76091-02	Medium Right	76091-08
Large Right	76091-04	Large Right	76091-10
Small Left	76091-01	Small Left	76091-07
Medium Left	76091-03	Medium Left	76091-09
Large Left	76091-05	Large Left	76091-11



JOBST® LOWER LIMB RANGES

SELECTING THE MOST APPROPRIATE JOBST® READY-TO-WEAR PRODUCT



JOBST UltraSheer Stockings

An exceptional combination of comfort, style and graduated compression therapy. JOBST® UltraSheer offers a fashionable sheer look in 4 styles and 3 colours to suit every patients needs.

- Maternity available
- Open toe available (only in knee high 15-20mmHg)

15 - 20 mmHg

20 - 30 mmHg

30 - 40 mmHg

Natural

Pattern

Black



JOBST forMen Casual Socks (unisex)

A soft, contemporary, relaxed design offering graduated compression. Dri release delivers superior moisture wicking technology to keep feet dry and cool. With its reinforced heel and strong yarns, JOBST® forMen Casual Socks are able to resist daily wear and tear. Contains a wool blend.

15 - 20 mmHg

20 - 30 mmHg

30 - 40 mmHg

Khaki

Black



JOBST forMen

With a unique knit structure, the ribbed design offers a fashionable men's look dress sock. JOBST forMen offers therapy and style with unparalleled comfort and durability.

15 - 20 mmHg

20 - 30 mmHg

30 - 40 mmHg

Black

JOBST® LOWER LIMB RANGES

SELECTING THE MOST APPROPRIATE JOBST® READY-TO-WEAR PRODUCT



JOBST ActiveWear Socks (unisex)

A soft, latex free, easy to apply medical compression sock that fits active lifestyles by providing the look and feel of athletic socks. The soft and seamless fit facilitates easy application and maximum comfort.

15 - 20 mmHg 20 - 30 mmHg 30 - 40 mmHg

White Black



JOBST Sport

The athletic compression socks for fast moisture transportation. Soft and lightweight with a breathable knit, they provide long wearing comfort for all activities.

15 - 20 mmHg

Pink/Grey Royal/Grey



JOBST Relief






























Graduated compression in a classical stocking appearance. With its remarkable durability and comfort, the product will offer lasting support and longevity.

- Open and closed toe available

15 - 20 mmHg 20 - 30 mmHg 30 - 40 mmHg

Beige Black

READY-TO-WEAR MEDICAL LEGWEAR

Name and Style	S		M		L		XL		Compression	Colour
	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe		
UltraSheer Stockings Knee High	75165-30	75165-25	75166-30	75166-25	75167-30	75167-25	75168-30	75168-25	15-20 mmHg	Natural Regular 
		75165-04		75166-04		75167-04		75168-04		Black 
		73633-00		73633-01		73633-02		73633-03		Diamond Pattern Black 
		75201-19		75202-19		75203-19		75204-19		Natural Regular 
		75201-04		75202-04		75203-04		75204-04		Black 
		73634-00		73634-01		73634-02		73634-03		Diamond Pattern Black 
UltraSheer Stockings Thigh High		75241-21		75242-21		75243-21		75244-21	30-40 mmHg	Natural 
		75241-04		75242-04		75243-04		75244-04		Black 
		75169-22		75170-22		75171-22		75172-22		Natural
		75169-04		75170-04		75171-04		75172-04		Black 
		73635-00		73635-01		73635-02		73635-03		Diamond Pattern Black 
		75205-23		75206-23		75207-23		75208-23		Natural 
UltraSheer SENSITIVE Stockings Thigh High		75205-04		75206-04		75207-04		75208-04	20-30 mmHg	Black 
		73636-00		73636-01		73636-02		73636-03		Diamond Pattern Black 
		75245-20		75246-20		75247-20		75248-20		Natural 
		75245-04		75246-04		75247-04		75248-04		Black 
		77686-12		77686-13		77686-14		77686-15		Natural 
		77686-04		77686-05		77686-06		77686-07		Black 
UltraSheer Stockings Waist High		7687-12		7687-13		7687-14		7687-15	15-20 mmHg	Natural 
		7687-04		7687-05		7687-06		7687-07		Black 
		77688-12		77688-13		77688-14		77688-15		Natural 
		77688-04		77688-05		77688-06		77688-07		Black 
		75173-12		75174-12		75175-12		75176-12		Natural 
		75173-04		75174-04		75175-04		75176-04		Black 
UltraSheer Stockings		73637-00		73637-01		73637-02		73637-03	20-30 mmHg	Diamond Pattern Black 
		75209-13		75210-13		75211-13		75212-13		Natural 
		75209-04		75210-04		75211-04		75212-04		Black 
		73638-00		73638-01		73638-02		73638-03		Diamond Pattern Black 
		75249-13		75250-13		75251-13		75252-13		Natural 
		75249-04		75250-04		75251-04		75252-04		Black 

DAILY

JOBST® SELECTION CHART

	75177-04	75178-04	75179-04	75180-04	75180-04	75180-01	15-20 mmHg	Natural
UltraSheer Stockings Maternity	75177-01	75178-01	75179-01	75180-01		75180-01	15-20 mmHg	Black
	75213-06	75214-06	75215-06				20-30 mmHg	Natural
	75213-03	75214-03	75215-03				20-30 mmHg	Black
for Men Casual Socks Knee High	75476-02	75477-02	75478-02			75479-02	15-20 mmHg	Black
	75476-03	75477-03	75478-03			75479-03	15-20 mmHg	Khaki (Beige)
	75484-02	75485-02	75486-02			75487-02	20-30 mmHg	Black
	75484-03	75485-03	75486-03			75487-03	20-30 mmHg	Khaki (Beige)
	75492-02	75493-02	75494-02			75495-02	30-40 mmHg	Black
	75492-03	75493-03	75494-03			75495-03	30-40 mmHg	Khaki (Beige)
	75253-01	75254-01	75255-01			75256-01	15-20 mmHg	Black
	75257-01	75258-01	75259-01			75260-01	20-30 mmHg	Black
for Men Socks Knee High	75265-01	75266-01	75267-01			75268-01	30-40 mmHg	Black
	75285-00	75286-00	75287-00			75288-00	15-20 mmHg	Black
	75261-00	75262-00	75263-00			75264-00	20-30 mmHg	Black
	75269-00	75270-00	75271-00			75272-00	30-40 mmHg	Black
	75147-03	75148-03	75149-03			75150-03	15-20 mmHg	Black
ActiveWear Socks Knee High	75147-02	75148-02	75149-02			75150-02	15-20 mmHg	White
	75155-03	75156-03	75157-03			75158-03	20-30 mmHg	Black
	75155-02	75156-02	75157-02			75158-02	20-30 mmHg	White
	75411-03	75412-03	75413-03			75414-03	30-40 mmHg	Black
	75411-02	75412-02	75413-02			75414-02	30-40 mmHg	White
	75289-80	75289-81	75289-82			75289-83	15-20 mmHg	Royal Blue
	75289-60	75289-61	75289-62			75289-63	15-20 mmHg	Pink
Relief Garments Knee High	75339-12	75340-12	75341-12			75342-12	15-20 mmHg	Black
	75339-11	75340-11	75341-11	75343-12	75343-11	75344-11	15-20 mmHg	Beige
	78050-04	78051-04	78052-04	75343-10	75343-10	75344-10	20-30 mmHg	Beige
	78054-04	78055-04	78056-04	77975-15	77975-16	77976-15	20-30 mmHg	Beige
	78042-08	78043-08	78044-08	77977-15	77977-16	77978-15	30-40 mmHg	Beige
	78059-08	78047-08	78048-08	78045-08	78045-09	78049-09	20-30 mmHg	Beige
	78476-10	78477-10	78478-10	78479-10	78479-11	78479-11	30-40 mmHg	Beige
	78480-10	78481-10	78482-10	78483-10	78483-10	78483-10	30-40 mmHg	Beige Left
	78031-05	78032-05	78033-05	78034-05	78034-05	78034-05	20-30 mmHg	Beige Right
	78031-06	78032-06	78033-06	78034-06	78034-06	78034-06	20-30 mmHg	Beige Left
Relief Garments Chaps	78035-05	78036-05	78037-05	78041-05	78041-05	78041-05	30-40 mmHg	Beige Left
	78035-06	78036-06	78037-06	78041-06	78041-06	78041-06	30-40 mmHg	Beige Right

WORK WEAR

SPORT

TRADITIONAL

JOBST® LOWER LIMB RANGES SIZING CHART

JOBST® ULTRASHEER, FORMEN, FORMEN CASUAL, ACTIVEWEAR, SPORT AND RELIEF SIZING



MEASUREMENTS TABLE

Style	Measurement	Small	Medium	Large	Large FC	X Large	X Large FC
Unisex	Ankle	18 - 21 cm	21 - 25 cm	25 - 29 cm	25 - 30 cm	29 - 33 cm	30 - 36 cm
	Calf	28 - 38 cm	30 - 42 cm	32 - 46 cm	46 - 61 cm	34 - 50 cm	46 - 61 cm
	Thigh	40 - 62 cm	46 - 70 cm	54 - 78 cm		60 - 81 cm	
	Hip	71 - 117 cm	76 - 127 cm	81 - 137 cm		102 - 166 cm	
Maternity	Ankle	18 - 21 cm	21 - 25 cm	25 - 29 cm		29 - 33 cm	
	Thigh	40 - 62 cm	46 - 70 cm	54 - 78 cm		60 - 81 cm	
	Hip	81 - 132 cm	86 - 142 cm	91 - 152 cm		112 - 180 cm	

JOBST® ULCERCARE™

JOBST® ULCERCARE™




JOBST® UlcerCARE™ enables effective management of hard to heal wounds on the lower leg while offering ease of use for patients and care givers - the perfect combination for providing patients with an improved quality of life.

JOBST® UlcerCARE™

Compression System:

- Provides graduated compression of 40 mmHg, the recommended pressure for the healing of these types of wounds
- Is a 2-in-1 stocking system designed for patient's self management to facilitate easy changing of the wound dressing and to allow for improved personal hygiene (e.g. bathing and skin care)
- JOBST UlcerCARE Liner provides mild compression (17mmHg) for the non-ambulatory patient to help manage lower leg oedema during bed rest
- Outer stocking comes with a zipper on the left or right side for comfortable closure on the opposite side of the wound. It is removed at night and applied first thing in the morning

MEASUREMENTS TABLE



Size	Ankle	Calf	Zipper on Left Side	Zipper on Right Side	JOBST UlcerCARE Replacement Liners (3 liners)
Small	18-21 cm	29-35 cm	73631-21	73631-61	73632-21
Medium	21-24 cm	34-41 cm	73631-22	73631-62	73632-22
Large	24-27 cm	40-47 cm	73631-23	73631-63	73632-23
X-Large	27-30 cm	46-54 cm	73631-24	73631-64	73632-24
XX-Large	30-33 cm	48-56 cm	73631-25	73631-65	73632-25
XXX-Large	33-36 cm	50-58 cm	73631-26	73631-66	73632-26
XXXX-Large	36-39 cm	52-60 cm	73631-27	73631-67	73632-27

JOBST® ELVAREX SOFT

JOBST® ELVAREX SOFT KNEE HIGH COMPRESSION GARMENTS




Ready to Wear Soft Flat Knit Compression Garments

JOBST® Elvarex® Soft Knee High garments are the perfect alternative for patients where management with circular knitted garments are no longer suitable.

They are indicated for the management of mild Lymphoedema and vascular insufficiency

- **Soft and Skin Friendly:** The soft knitted yarn combines skin protecting softness with excellent breathability, making it an ideal solution for those with sensitive and fragile skin
- **Therapeutically effective:** Combining the right level of compression with the right fabric stiffness, JOBST Elvarex Soft offers an extremely comfortable and effective treatment for low to moderate oedema



MEASUREMENTS TABLE

JOBST ELVAREX SOFT 20-30 MMHG SILICONE BAND - SHORT 2.5 CM, LONG 5 CM

Size	I	I plus	II	II plus	III	III plus
Long (ID AD = 40-45)	79559-012	79559-018	79559-013	79559-019	79559-014	79559-020
Short (ID AD = 34-39)	79558-012	79558-018	79558-013	79558-019	79558-014	79558-020
Size	IV	IV plus	V	V plus	VI	VI plus
Long (ID AD = 40-45)	79559-015	79559-021	79559-016	79559-022	79559-017	79559-023
Short (ID AD = 34-39)	79558-015	79558-021	79558-016	79558-022	79558-017	79558-023

JOBST ELVAREX SOFT 30-40 MMHG SILICONE BAND - SHORT 2.5 CM, LONG 5 CM

Size	I	I plus	II	II plus	III	III plus
Long (ID AD = 40-45)	73552-012	73552-018	73552-013	73552-019	73552-014	73552-020
Short (ID AD = 34-39)	73579-012	73579-018	73579-013	73579-019	73579-014	73579-020
Size	IV	IV plus	V	V plus	VI	VI plus
Long (ID AD = 40-45)	73552-015	73552-021	73552-016	73552-022	73552-017	73552-023
Short (ID AD = 34-39)	73579-015	73579-021	73579-016	73579-022	73579-017	73579-023
cD	28-32	30-34	32-36	34-38	36-40	38-42
cD plus	32-36	34-38	36-40	38-42	40-44	42-46
cC	29-33	31-35	33-37	35-39	37-41	39-43
cC plus	33-37	35-39	37-41	39-43	41-45	43-47
cB1	23-26	25-28	27-30	29-32	31-34	33-36
cB1 plus	24-27	26-29	28-31	30-33	32-35	34-37
cB	18-20	20-22	22-24	24-26	26-28	28-30
cY	26-31	28-33	29-35	31-37	32-38	33-40
cA	17-22	19-24	21-26	23-29	25-32	27-34

THE JOBST® FarrowWrap® PRODUCT RANGE

LITE

Available variants

- Upper and lower limb
- Ready-to-wear (RTW)
- Sleeve • Gauntlet
- Thigh • Leg
- Foot • Knee



A double laminated fabric suitable for mild to moderate compression needs, making it especially suitable for palliative and paediatric patients. The fabric is both durable and breathable.

- Compression: 20-30 mmHg
- Mild to moderate chronic oedema/lymphoedema
- Conforms to limb shape
- Good choice for patients who have neuropathy or peripheral vascular disease (PVD)
- Recommended for paediatric patients
- Suitable for sensitive skin and as such palliative care patients
- Machine washable and can be tumble dried on low heat

STRONG

Available variants

- Lower limb
- Ready-to-wear (RTW)
- Thigh
- Leg
- Foot
- Knee



Premium fabric with strong compression and soft inner layer suitable for patients with a good limb shape and moderate to severe oedema.

- Compression: 30-40 mmHg
- Moderate to severe chronic oedema/lymphoedema with a good limb shape
- Less rigid than CLASSIC
- Contoured bands conform to the limb and resist sliding
- Highly breathable and durable
- Suitable for minor skin folds and distortion
- Material less likely to pill or fray
- Machine washable and can be tumble dried on low heat

CLASSIC

Available variants

- Lower limb
- Ready-to-wear (RTW)
- Thigh
- Leg
- Foot



The original material that has been clinically proven over many years of use on thousands of patients worldwide. The rigid nature of the material offers performance beyond the compression range making it excellent for patients with moderate to severe compression levels.

- Compression: 30-40 mmHg
- Moderate to severe chronic oedema/lymphoedema
- Suitable for oedema with skin folds and irregular limb shapes
- Durable, rigid fabric with a comfortable feel
- Reinforced spine for easier donning
- Less likely to roll down and will not work into skin folds
- Hand wash only and lay flat to dry

THE JOBST® FarrowWrap®

PRODUCT ASSORTMENT

Mild to Moderate 20-30mmHg


JOBST® FarrowWrap® LITE

JOBST® FarrowHybrid® ADI

Moderate to Strong 30-40mmHg

JOBST® FarrowWrap® STRONG

JOBST® FarrowWrap® CLASSIC



Footpiece	Legpiece	Thighpiece	FarrowHybrid Foot Liner	Hand Gauntlet Foam Squares	Toe Cap/Glove
LITE RTW	LITE RTW	LITE RTW	ADI RTW	LITE RTW	LITE RTW
STRONG RTW	STRONG RTW	STRONG RTW			
CLASSIC RTW	CLASSIC RTW	CLASSIC RTW			

ACCESSORIES

JOBST® Farrow GarmentGrip™ Roll



JOBST® Farrow GarmentGrip can be used over the top of the wrap/garment with part of it against the skin and the other part overlapping the garment to help prevent slippage. Garment Grip can be cut to size to fit the shape and contour of the body part and each roll comes with several VELCRO® Brand fasteners to secure it in place. Jobst Farrow Garment Grip can also be used to add extra compression over a wrap/garment where needed. It is washable and therefore reusable. It is available in differing roll lengths and extra wide widths to suit all limb parts.

Code	Size	Limb	Colour
76672-05	10cm x 60cm	arm	Beige
76672-00	14cm x 80cm	thigh	Beige
76672-01	14cm x 120cm	thigh long	Beige
76672-03	20cm x 120cm	wide	Beige
76672-04	30cm x 150cm	extra wide	Beige
76672-02	14cm x 550cm	extra long	Beige


THE JOBST® FarrowWrap® HOW TO MEASURE

RTW FOOTPIECE SIZE CHART

READY-TO-WEAR FOOTPIECE

		Xsmall	Small	Medium	Large	X Large
A'	Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
X	Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm
X	Long	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm

The length is measured from the most posterior aspect of the heel to the base of the 5th toe.




RTW LEGPIECE SIZE CHART

READY-TO-WEAR LEGPIECE

		Xsmall	Small	Medium	Large	X Large
C	Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
B	Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D	Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm
A-D	Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm


Length A-D = Floor to 2 finger widths below crease of knee. Measure posteriorly FOLLOWING contour of leg.



RTW THIGHPIECE SIZE CHART

Size	G Groin	E Above Knee	E-G Short	E-G Regular	E-G Tall
XS	60-70 cm	47-52 cm	16-20 cm	21-26 cm	27-31 cm
S	65-75 cm	53-59 cm	16-20 cm	21-26 cm	27-31 cm
M	70-80 cm	60-66 cm	16-20 cm	21-26 cm	27-31 cm
L	75-90 cm	67-73 cm	16-20 cm	21-26 cm	27-31 cm
XL	85-100 cm	74-80 cm	16-20 cm	21-26 cm	27-31 cm


Length E-G = Measure **MEDIALLY** from just above patella to the groin. Follow outer contour of any skin folds/lobules.



THE JOBST® FARROW LOWER EXTREMITY


ORDERING INFORMATION

JOBST® FARROWWRAP® LITE (20-30MMHG) FOOTPIECE




Size	Product Code	Length	Colour
XS	76652-05	REGULAR	TAN
S	76652-06	REGULAR	TAN
M	76652-07	REGULAR	TAN
L	76652-08	REGULAR	TAN
XL	76652-09	REGULAR	TAN
XS	76652-00	LONG	TAN
S	76652-01	LONG	TAN
M	76652-02	LONG	TAN
L	76652-03	LONG	TAN
XL	76652-04	LONG	TAN

JOBST® FARROWWRAP® LITE (20-30MMHG) LEGPIECE



Size	Product Code	Length	Colour
XS	76654-00	REGULAR	TAN
S	76654-01	REGULAR	TAN
M	76654-02	REGULAR	TAN
L	76654-03	REGULAR	TAN
XL	76654-04	REGULAR	TAN
XS	76654-05	TALL	TAN
S	76654-06	TALL	TAN
M	76654-07	TALL	TAN
L	76654-08	TALL	TAN
XL	76654-09	TALL	TAN

JOBST® FARROWWRAP® LITE (20-30MMHG) THIGHPIECE




Size	Product Code	Length	Colour
XS	76655-05	SHORT	TAN
S	76655-06	SHORT	TAN
M	76655-07	SHORT	TAN
L	76655-08	SHORT	TAN
XL	76655-09	SHORT	TAN
XS	76655-00	REGULAR	TAN
S	76655-01	REGULAR	TAN
M	76655-02	REGULAR	TAN
L	76655-03	REGULAR	TAN
XL	76655-04	REGULAR	TAN
XS	76655-10	TALL	TAN
S	76655-11	TALL	TAN
M	76655-12	TALL	TAN
L	76655-13	TALL	TAN
XL	76655-14	TALL	TAN

The JOBST® FARROW LOWER EXTREMITY


ORDERING INFORMATION

JOBST® FARROWWRAP® STRONG (30-40MMHG) FOOTPIECE




Size	Product Code	Length	Colour
XS	76646-04	REGULAR	TAN
S	76646-05	REGULAR	TAN
M	76646-06	REGULAR	TAN
L	76646-07	REGULAR	TAN
XL	76646-08	REGULAR	TAN
XS	76646-09	LONG	TAN
S	76646-00	LONG	TAN
M	76646-01	LONG	TAN
L	76646-02	LONG	TAN
XL	76646-03	LONG	TAN

JOBST® FARROWWRAP® STRONG (30-40MMHG) LEGPIECE



Size	Product Code	Length	Colour
XS	76648-00	REGULAR	TAN
S	76648-01	REGULAR	TAN
M	76648-02	REGULAR	TAN
L	76648-03	REGULAR	TAN
XL	76648-04	REGULAR	TAN
XS	76648-05	TALL	TAN
S	76648-06	TALL	TAN
M	76648-07	TALL	TAN
L	76648-08	TALL	TAN
XL	76648-09	TALL	TAN

JOBST® FARROWWRAP® STRONG (30-40MMHG) THIGHPIECE




Size	Product Code	Length	Colour
XS	76649-05	SHORT	TAN
S	76649-06	SHORT	TAN
M	76649-07	SHORT	TAN
L	76649-08	SHORT	TAN
XL	76649-09	SHORT	TAN
XS	76649-00	REGULAR	TAN
S	76649-01	REGULAR	TAN
M	76649-02	REGULAR	TAN
L	76649-03	REGULAR	TAN
XL	76649-04	REGULAR	TAN
XS	76649-10	TALL	TAN
S	76649-11	TALL	TAN
M	76649-12	TALL	TAN
L	76649-13	TALL	TAN
XL	76649-14	TALL	TAN

The JOBST® FARROW LOWER EXTREMITY


ORDERING INFORMATION

JOBST® FARROWWRAP® CLASSIC (30-40MMHG) FOOTPIECE




Size	Product Code	Length	Colour
XS	76656-05	REGULAR	TAN
S	76656-06	REGULAR	TAN
M	76656-07	REGULAR	TAN
L	76656-08	REGULAR	TAN
XL	76656-09	REGULAR	TAN
XS	76656-00	LONG	TAN
S	76656-01	LONG	TAN
M	76656-02	LONG	TAN
L	76656-03	LONG	TAN
XL	76656-04	LONG	TAN

JOBST® FARROWWRAP® CLASSIC (30-40MMHG) LEGPIECE



Size	Product Code	Length	Colour
XS	76657-00	REGULAR	TAN
S	76657-01	REGULAR	TAN
M	76657-02	REGULAR	TAN
L	76657-03	REGULAR	TAN
XL	76657-04	REGULAR	TAN
XS	76657-05	TALL	TAN
S	76657-06	TALL	TAN
M	76657-07	TALL	TAN
L	76657-08	TALL	TAN
XL	76657-09	TALL	TAN

JOBST® FARROWWRAP® CLASSIC (30-40MMHG) THIGHPIECE



Size	Product Code	Length	Colour
XS	76658-06	SHORT	TAN
S	76658-07	SHORT	TAN
M	76658-08	SHORT	TAN
L	76658-09	SHORT	TAN
XL	76658-10	SHORT	TAN
XS	76658-00	REGULAR	TAN
S	76658-01	REGULAR	TAN
M	76658-02	REGULAR	TAN
L	76658-03	REGULAR	TAN
XL	76658-04	REGULAR	TAN
XS	76658-12	TALL	TAN
S	76658-13	TALL	TAN
M	76658-14	TALL	TAN
L	76658-15	TALL	TAN
XL	76658-16	TALL	TAN

THE JOBST® FarrowWrap® LOWER EXTREMITY

KNEEPIECE




Adjustable knee piece, easy to apply & comfortable to wear.

Can be worn in combination with JOBST® FarrowWrap® Legpieces and Thighpieces, or as a single unit (e.g. in combination with an Elvarex knee high stocking)

Available in two fabrics: **LITE** and **STRONG**

JOBST® FARROWWRAP® KNEEPIECE SIZING CHART

	Xsmall	Small	Medium	Large	X Large
5cm above patella	47 - 52 cm	53 - 59 cm	60 - 66 cm	67 - 73 cm	74 - 80 cm
Calf	36 - 43 cm	42 - 50 cm	48 - 58 cm	53 - 63 cm	58 - 68 cm




ORDERING INFORMATION

JOBST® FARROWWRAP® LITE (20-30MMHG) KNEEPIECE

Size	Product Code	Colour
XS	73367-00	TAN
S	73367-01	TAN
M	73367-02	TAN
L	73367-03	TAN
XL	73367-04	TAN

JOBST® FARROWWRAP® STRONG (30-40MMHG) KNEEPIECE

Size	Product Code	Colour
XS	73368-00	TAN
S	73368-01	TAN
M	73368-02	TAN
L	73368-03	TAN
XL	73368-04	TAN



THE JOBST® FARROW LOWER EXTREMITY

TOE CAP




Ready to wear glove which can be worn on either right or left foot.

Fine seams, with ultra thin fabric, allows air to pass to enhance patient comfort. The toes can be easily trimmed with scissors to exactly the required length without unravelling.

MEASUREMENTS TABLE

Size	cA
XS	20-22 cm
S	22-24 cm
M	24-26 cm
L	26-28 cm



ORDERING INFORMATION

Size	Product Code	Compression	Colour
XS	72947-00	20 - 30 mmHg	Beige
S	72947-01	20 - 30 mmHg	Beige
M	72947-02	20 - 30 mmHg	Beige
L	72947-03	20 - 30 mmHg	Beige



JOBST® FarrowHybrid® LINERS

LINERS & ACCESSORIES



FarrowHybrid® is a patented foot compression garment that provides compression in the foot up to the ankle.

- From the ankle and above, it serves as a liner
- Available in standard and wide sizes
- Great alternative to a footpiece
- Low profile, easily fits into a shoe

	Male Shoe Size	Female Shoe Size
S	3 - 6	4 - 7
M	7 - 9	8 - 10
L	10 - 12	11 - 13
XL	13 - 15	14 - 16

Hybrid Liners

Measure the mid foot A circumference. Two mid foot sizes available: 20-27 cm (regular foot) and 28-35 cm (wide foot). Request the patient's shoe size (US Sizing).

JOBST® FARROWHYBRID KNEE HIGH LINERS (20-30MMHG)

Size	Code	Fabric	Female shoe size	Male shoe size	Colour	Mid foot Circumference
S	76664-00	Hybrid ADI	US Female Shoe 4-7	US Male Shoe 3-6	Beige	Mid foot 20-27 cm
M	76664-01	Hybrid ADI	US Female Shoe 8-10	US Male Shoe 7-9	Beige	Mid foot 20-27 cm
L	76664-02	Hybrid ADI	US Female Shoe 11-13	US Male Shoe 10-12	Beige	Mid foot 20-27 cm
XL	76664-03	Hybrid ADI	US Female Shoe 14-16	US Male Shoe 13-15	Beige	Mid foot 20-27 cm
Size	Code	Fabric	Female shoe size	Male shoe size	Colour	Mid foot Circumference
S	76664-12	Hybrid ADI	US Female Shoe 4-7	US Male Shoe 3-6	Beige	Mid foot 28-35 cm
M	76664-13	Hybrid ADI	US Female Shoe 8-10	US Male Shoe 7-9	Beige	Mid foot 28-35 cm
L	76664-14	Hybrid ADI	US Female Shoe 11-13	US Male Shoe 10-12	Beige	Mid foot 28-35 cm
XL	76664-15	Hybrid ADI	US Female Shoe 14-16	US Male Shoe 13-15	Beige	Mid foot 28-35 cm
Size	Code	Fabric	Female shoe size	Male shoe size	Colour	Mid foot Circumference
S	76664-16	Hybrid ADI	US Female Shoe 4-7	US Male Shoe 3-6	Black	Mid foot 28-35 cm
M	76664-17	Hybrid ADI	US Female Shoe 8-10	US Male Shoe 7-9	Black	Mid foot 28-35 cm
L	76664-18	Hybrid ADI	US Female Shoe 11-13	US Male Shoe 10-12	Black	Mid foot 28-35 cm
XL	76664-19	Hybrid ADI	US Female Shoe 14-16	US Male Shoe 13-15	Black	Mid foot 28-35 cm

Lower Limb Liners

JOBST® KNEE HIGH LINERS - SILVER CLOSED TOE

Size	Code	Style	Calf Circumference	Colour
Standard	76666-21	Closed Toe	Up to 58cm calf	Beige
Wide	76666-22	Closed Toe	52 - 70cm calf	Beige

JOBST® KNEE HIGH LINERS - TERRY CLOTH

Size	Code	Style	Calf Circumference	Colour
S	76669-00	Open Toe	< 40 cm calf	Beige
M	76669-02	Open Toe	40 - 70cm calf	Beige
L	76669-04	Open Toe	70 - 120cm calf	Beige

JOBST® THIGH HIGH LINERS - TERRY CLOTH

Size	Code	Style	Calf Circumference	Colour
S	76670-00	Open Toe	< 40 cm calf	Beige
M	76670-02	Open Toe	40 - 70cm calf	Beige
L	76670-04	Open Toe	70 - 120cm calf	Beige

JOBST FarrowWrap® UPPER EXTREMITY

JOBST® FARROWWRAP® ARMPIECE - SUITABLE FOR 24-HOUR WEAR



LITE RTW Armpiece

LITE RTW Armpiece is an anatomically contoured effective solution to control upper extremity oedema. It is easy to apply and remove with one hand and designed for mild to moderate oedema control. Light and comfortable with just one layer but with a robust short-stretch performance.

- 20-30 mmHg compression

ORDERING INFORMATION

JOBST® FARROWWRAP® LITE (20-30MMHG) ARMPIECE LONG


Option/Colour	Size	Code	Option/Colour	Size	Code
Left/Tan	S	76650-35	Right/Tan	S	76650-00
Left/Tan	M	76650-01	Right/Tan	M	76650-02
Left/Tan	L	76650-03	Right/Tan	L	76650-04
Left/Black	S	76650-17	Right/Black	S	76650-18
Left/Black	M	76650-19	Right/Black	M	76650-20
Left/Black	L	76650-21	Right/Black	L	76650-22

JOBST® FARROWWRAP® LITE (20-30MMHG) ARMPIECE REGULAR

Option/Colour	Size	Code	Option/Colour	Size	Code
Left/Tan	S	76650-05	Right/Tan	S	76650-06
Left/Tan	M	76650-07	Right/Tan	M	76650-08
Left/Tan	L	76650-09	Right/Tan	L	76650-10
Left/Black	S	76650-23	Right/Black	S	76650-24
Left/Black	M	76650-25	Right/Black	M	76650-26
Left/Black	L	76650-27	Right/Black	L	76650-28

JOBST® FARROWWRAP® LITE (20-30MMHG) ARMPIECE SHORT

Option/Colour	Size	Code	Option/Colour	Size	Code
Left/Tan	S	76650-11	Right/Tan	S	76650-12
Left/Tan	M	76650-13	Right/Tan	M	76650-14
Left/Tan	L	76650-15	Right/Tan	L	76650-16
Left/Black	S	76650-29	Right/Black	S	76650-30
Left/Black	M	76650-31	Right/Black	M	76650-32
Left/Black	L	76650-33	Right/Black	L	76650-34



Short	40-43 cm
Regular	44-47 cm
Long	48-51 cm

NOTE: Measure length along the lateral outer arm from the wrist to the axilla with a slight bend in the arm.

SIZE CHART

Wrist	S			M			L		
	14-18 cm	16-21 cm	19-25 cm	20-27 cm	25-34 cm	30-40 cm			
Elbow	S			M			L		
	20-27 cm			25-34 cm			30-40 cm		
Axilla	S			M			L		
	22-31 cm			29-39 cm			32-45 cm		

JOBST FarrowWrap® UPPER EXTREMITY

JOBST® FARROWWRAP® HAND - SUITABLE FOR 24-HOUR WEAR




JOBST FarrowWrap LITE Hand Gauntlet 20-30 mmHg

JOBST FarrowWrap LITE Hand Gauntlet is ambidextrous and comes with 4 mm and 8 mm foam squares for added compression to the dorsal hand. Suitable for mild to severe oedemas.

Benefits that give your patient greater freedom to enjoy life:

- Foam inserts help reduce oedema that can pool on the top of the hand
- Ambidextrous, i.e. can be used on either hand
- Durable and comfortable

SIZE CHART



JOBST® FARROWWRAP® LITE (20-30MMHG) HAND GAUNTLET

	A	B	C
XS	16-18 cm	17-19 cm	15-18 cm
S	18-20 cm	19-21 cm	17-20 cm
M	20-22 cm	21-23 cm	19-22 cm
L	22-24 cm	23-26 cm	22-24 cm

JOBST® FARROWWRAP LITE (20-30MMHG) HAND GAUNTLET

Size	Code	Fabric	Style	Limb	Colour
XS	76653-00	Lite	Hand (with foam)	Ambidextrous	Beige
S	76653-01	Lite	Hand (with foam)	Ambidextrous	Beige
M	76653-02	Lite	Hand (with foam)	Ambidextrous	Beige
L	76653-03	Lite	Hand (with foam)	Ambidextrous	Beige

Clinical Equipment

We can supply you with a wide range of clinical equipment. From treatment tables to massage chairs, trolleys to stools and everything in between.

Visit opchealth.com.au to view our entire range



BELLISSE® BRA AND PADS

BECAUSE A HUG SHOULD NEVER HURT




Bellisse®

Bellisse® is a medical compression garment that facilitates drainage of excess lymphatic fluid in the chest and breast area following breast cancer, and other chest or truncal surgeries.

This garment will provide you with:

- Gentle compression of the torso while directing the excess lymphatic load to alternate functional pathways
- The use of specialty pads is simplified by addressing specific areas of breast oedema or chest wall fibrosis
- Available in a large range of sizes, with unprecedented adjustability: a front zipper for ease of donning, and a variable position back closure - allowing for optimal fit
- The unique, non-restrictive design is comfortable, day and night, and aids in pain reduction and swelling following surgery.

Colour	Fabric
Buff	Polyester/Nylon/Spandex



When do you need to add a breast pad?

Breast pads are important when there is pitting oedema (when finger pressure leaves a lingering indentation in the skin) and/or fibrosis (stiffness or scarring) resulting from surgery or radiation. Both conditions impede the natural flow of fluid and contribute to chronic swelling in the breast or chest wall.

Axilla Pad

The Axilla Pad was created in response to a request by a breast surgeon looking for a "tapeless solution" immediately following axillary node dissection. The flat pad is placed against the chest wall and held in place by the Bellisse®.




Colour	Fabric
Buff	Organic Cotton & Spandex with elastic band
Code	Size
75925-02	Axila Pad 30-38 cm Small/ Medium Buff
75925-03	Axila Pad 38-46 cm Large/ Xlarge Buff

The size is determined by the widest biceps measurement:
30-38 cm or 38-46 cm

BELLISSE® BRA AND PADS

BECAUSE A HUG SHOULD NEVER HURT



Lateral Pad

The Lateral Pad fits into the inside pocket of the Bellisse® and covers the Bellisse® or other compression garment. It is ideal for softening fibrotic areas on the lateral breast.

Colour	Fabric
Buff	Polartec® Power Dry®
Code	Size
76166-03	Lateral Pad Small Buff
76166-04	Lateral Pad Medium Buff
76166-05	Lateral Pad Large Buff

Post-Lumpectomy Pad

The Post-Lumpectomy Pad features a posterior extension that addresses swelling and/or radiation fibrosis over the affected breast and serratus anterior. This pad covers the affected breast completely and extends to the spine.

Colour	Fabric
Buff	Polartec® Power Dry®
Code	Size
75931-04	Post Lumpectomy Pad A/B Cup Small Buff
75931-05	Post Lumpectomy Pad C Cup Medium Buff
75931-06	Post Lumpectomy Pad D/DD Cup Large Buff
75931-07	Post Lumpectomy Pad E Cup Xlarge Buff

BELLISSE® BRA SIZE CHART

HOW TO MEASURE

To determine your patient's band size:


- 1a. Measure circumferentially below the bust where a regular bra finishes.
- 1b. Drop the tape 5 - 7 cm below this point and measure around the bottom of the rib cage.

Use the larger of the two rib cage measurements for your patient's band size. Round odd numbers up.

To determine your patient's cup size:

2. Measure circumferentially around the fullest part of the bust.
3. Subtract largest ribcage measurement from the bust measurement to determine cup size.
 - a. If the difference is 5 cm or less between 1 and 2, select cup size A/B
 - b. If the difference is 5 - 10 cm between 1 and 2, select cup size C/D
 - c. If the difference is 10 - 15 cm between 1 and 2, select cup size DD/E
 - d. If the difference is more than 15 cm, then go up one band size

Size	Centimeters
30	76.3
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127



Difference	Cup Size
5cm or less	A/B
5 - 10cm	C/D
10 - 15cm	DD/E

If more than 15cm, go up one band size. Bellisse Size:

JOVIPAK HEAD AND NECK COMPRESSION GARMENTS

JOVIPAK HEAD AND NECK COMPRESSION GARMENTS

Cancer diagnosis of the head and neck, which results in lymphoedema, is a lifelong management challenge for patients and therapists. Skilled treatment by qualified professionals is required to successfully decrease symptoms and prevent complications. JoViPak has specially designed a line of compression garments for this condition.

These garments are well accepted by therapists and are highly tolerated by patients.

The unique feature of this line of garments is the highly resilient shredded foam which creates high- and low-pressure points in the skin and underlying tissue. Continuous flow channeling guides excess fluid away from congested areas and into functional lymphatic transport vessels.

Every garment is adjustable, washable, re-usable, and easy to apply. All garments are Ready Made.



Chin Strap - Standard

The Chin Strap - Standard is designed to address oedema of the neck, mandible and the anterior/lateral aspects of the face (partially covers the ears). The open neck seam accommodates a choice of padding. The Velcro® straps at the neck and top of head can be adjusted to patient's comfort level, making the garment suitable for overnight wear

Colour	Fabric
Buff	Organic Cotton/Lycra®

Note: this garment does not have any padding. JOBST® JoViPak® recommends a Neck Pad - Universal with this product for maximum effectiveness.



Neck Pad - Universal

The Neck Pad - Universal is designed to be used with the Chin Strap - Standard, addressing bilateral head & neck oedema where there is swelling of the neck in front of the ear(s).

Colour	Fabric
Ivory	Organic Cotton & Spandex
Code	Size
76133-00	Chin Strap Standard Small Size Ivory
76133-01	Chin Strap Standard Medium Size Ivory
76133-02	Chin Strap Standard Large Size Ivory
76133-03	Chin Strap Standard Xlarge Size Ivory
76133-04	Chin Strap Standard XXlarge Size Ivory
77610-00	Neck Pad Universal Shape Small/ medium Ivory
77610-01	Neck Pad Universal Shape Large/ Xlarge Ivory

HEAD AND NECK SIZE CHART

Size	Head circumference	Neck Circumference	Chin to crown circumference
S	53.3 - 55.9	28 - 33	60 - 63.5
M	55.9 - 58.4	33 - 38.1	63.5 - 67.5
L	58.4 - 61	38.1 - 43.2	67.5 - 71
XL	61 - 63.5	43.2 - 48.3	71 - 75
2XL	63.5 - 66	48.3 - 52.1	75 - 78.5

Comprilan®

COMPRESSION BANDAGES



100% Cotton Short-Stretch Bandage

Comprilan® can be used for phlebological and lymphological indications to support the muscle pump action and improve the haemodynamics of venous return. This short-stretch bandage is durable and extremely resistant.

- Washable and re-usable
- Air permeable
- Bandage ends do not fray

Size	Code	Pack Qty	Case Qty
6 cm x 5 m	01026-00	1 x 5	40 x 5
8 cm x 5 m	01027-00	1 x 5	34 x 5
10 cm x 5 m	01028-00	1 x 5	28 x 5
12 cm x 5 m	01029-00	1 x 5	24 x 5

JOBST® Compriz

TWO-LAYER BANDAGE SYSTEM

JOBST® Compriz has accessible guidance tools to aid application and promote the correct compression. Two compression levels, ensuring the right level can be selected for both venous and mixed aetiology leg ulcers.

Two choices to achieve the compression needed. Latex Free.



Type of disease	Incidence	ABPI	Compression level	
Venous	~ 60%	>0.8	40 mmHg	JOBST Compriz
Mixed aetiology	~ 20%	0.5-0.8	20-30 mmHg	JOBST Compriz lite
Arterial	~ 20%	<0.5	No compression	

Code	Ankle Size	Compression Level	Quantity	
76271-00	18-25 cm	Regular 40 mmHg	1 Kit	Latex free
76271-01	25-32 cm	Regular 40 mmHg	1 Kit	Latex free
76271-02	18-25 cm	Lite 30 mmHg	1 Kit	Latex free
76271-03	25-32 cm	Lite 30 mmHg	1 Kit	Latex free

JOBST® Accessories

JOBST ROLL ON ADHESIVE



Jobst Roll On Adhesive

- is a roll-on body fixative used to hold JOBST® LegWear in place
- is pliable, moves with your body and can be washed off easily with water
- does not stain
- Available in a 59ml bottle

Code	79274-01
------	----------

JOBST DONNING GLOVES



Jobst Donning Gloves

Jobst Donning Gloves are designed to assist in donning or doffing compression stockings. These cotton knit gloves with latex coating help grip stockings without tearing or ripping them. Donning gloves also help protect the stockings and socks from punctures and runs that may be caused by finger nails or jewellery. Not recommended for sheer garments

Code	Size	Unit of Measure
79799-03	S	EA
79799-04	M	EA

JOBST STOCKING DONNER



Jobst Stocking Donner

Due to the nature of compression stockings, people may find putting them on to be initially a challenging task. The JOBST Stocking Donner, which features a lightweight frame and foam grip handles, is designed to aid in the application of all JOBST Medical LegWear.

Code	79668-00
------	----------



CORE-SPUN[®]


by THERAFIRM[®]



Ultra stretchy, breathable yarns. Soft, comfortable band is non-restrictive while staying in place all day.



Delivers lab tested true gradient compression to promote better blood flow, prevent mild to moderate swelling, and relieve tired, achy legs and feet.



Moisture wicking CoolMax[®]* fibers provide dry comfortable coolness.
*Coolmax is a trademark of Invista



Reinforced toe and heel for durability.



Smooth toe seam.



WHERE TO MEASURE

1. Measure the largest portion of the calf
2. Measure just above the ankle bone

KNIT-RITE CORE-SPUN

CORE-SPUN BY THERAFIRM®

Core-Spun by Therafirm Support Socks, suitable for both men and women, look and feel like a soft, comfortable everyday sock with the added benefit of true medical graded graduated compression.

The ultra-stretchy Core-Spun yarns make the sock easier to put on (don), whilst the sock thickness makes them suitable for all-day wear. The added benefit of utilising Core-Spun Yarns are the increased leg circumferences that can be accommodated.

The true medical grade graduated compression promotes better circulation and helps prevent oedema and leg fatigue.

Core-Spun by Therafirm are produced using latest technology moisture wicking yarns to promote moisture to be wicked away from the skin to the exterior of the sock, up and away from the foot area, therefore resulting in a sock that is comfortable and cool.

Available in Ready-To-Wear sizes Sml - XX/Large, Ccl. I or Ccl.II, short or regular lengths in either Khaki or Black (with X-Static) as standard, with other colours available upon request.

SIZE CHART

Size	Ankle	Calf	Length
S Short	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	16"(41cm) or less
S Regular	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	greater than 16"(41cm)
M Short	8"-10" (21-26cm)	12"-17.5" (31-45cm)	16"(41cm) or less
M Regular	8"-10" (21-26cm)	12"-17.5" (31-45cm)	greater than 16"(41cm)
L Short	9"-11.5" (23-29cm)	13"-19" (33-49cm)	16"(41cm) or less
L Regular	9"-11.5" (23-29cm)	13"-19" (33-49cm)	greater than 16"(41cm)
XL Short	11"-15" (28-38cm)	17"-23" (43-58cm)	16"(41cm) or less
XL Regular	11"-15" (28-38cm)	17"-23" (43-58cm)	greater than 16"(41cm)
XXL Short	13"-17" (33-43cm)	18.5"-25.5" (47-65cm)	16"(41cm) or less
XXL Regular	13"-17" (33-43cm)	18.5"-25.5" (47-65cm)	greater than 16"(41cm)

SHORT LENGTH

Size	Code	Compression	Colour
S	17511	15-20mmHg	Khaki
M	17521	15-20mmHg	Khaki
L	17531	15-20mmHg	Khaki
XL	17541	15-20mmHg	Khaki
XXL	17551	15-20mmHg	Khaki
S	17711	15-20mmHg	Black
M	17721	15-20mmHg	Black
L	17731	15-20mmHg	Black
XL	17741	15-20mmHg	Black
XXL	17751	15-20mmHg	Black
S	18511	20-30mmHg	Khaki
M	18521	20-30mmHg	Khaki
L	18531	20-30mmHg	Khaki
XL	18541	20-30mmHg	Khaki
XXL	18551	20-30mmHg	Khaki
S	18711	20-30mmHg	Black
M	18721	20-30mmHg	Black
L	18731	20-30mmHg	Black
XL	18741	20-30mmHg	Black
XXL	18751	20-30mmHg	Black

REGULAR LENGTH

Size	Code	Compression	Colour
S	17512	15-20mmHg	Khaki
M	17522	15-20mmHg	Khaki
L	17532	15-20mmHg	Khaki
XL	17542	15-20mmHg	Khaki
XXL	17552	15-20mmHg	Khaki
S	17712	15-20mmHg	Black
M	17722	15-20mmHg	Black
L	17732	15-20mmHg	Black
XL	17742	15-20mmHg	Black
XXL	17752	15-20mmHg	Black
S	18512	20-30mmHg	Khaki
M	18521	20-30mmHg	Khaki
L	18532	20-30mmHg	Khaki
XL	18542	20-30mmHg	Khaki
XXL	18552	20-30mmHg	Khaki
S	18712	20-30mmHg	Black
M	18722	20-30mmHg	Black
L	18732	20-30mmHg	Black
XL	18742	20-30mmHg	Black
XXL	18752	20-30mmHg	Black


THERAFIRM®



Ultra stretchy, breathable yarns. Soft, comfortable band is non-restrictive while staying in place all day.



Delivers lab tested true gradient compression to promote better blood flow, prevent mild to moderate swelling, and relieve tired, achy legs and feet.



Moisture wicking CoolMax®* fibers provide dry comfortable coolness.
*Coolmax is a trademark of Invista



Reinforced toe and heel for durability.



Smooth toe seam.



WHERE TO MEASURE

1. Measure the largest portion of the calf
2. Measure just above the ankle bone

KNIT-RITE THERAFIRM

MEN'S TROUSER SOCK

Stylish below knee dress sock featuring ribbed style true graduated compression featuring non-binding comfort top to keep the sock in place all day. Features a comfort heel & toe. Sheer two way stretch nylon & spandex. Designed to prevent oedema, leg discomfort and deep vein thrombosis for long distant travellers.

Circular knitted, available in Ready-To-Wear sizes Sml - X/Large, Ccl.I or Ccl.II, in either Khaki/Black as standard, with other colours available upon special request.

SIZE CHART

Size	Ankle	Calf	Length	Size	Ankle	Calf	Length
S	7"-9" (18-23cm)	10.50"-15" (27-38cm)	11"-14" (28-36cm)	L	10"-14" (25-36cm)	12"-17" (31-43cm)	13"-16" (33-41cm)
M	9"-11" (18-23cm)	11"-16.50" (28-42cm)	12"-15" (31-38cm)	XL	11"-15" (28-38cm)	13"-19" (33-48cm)	14"-18" (36-46cm)

Size	Code	Compression	Colour	Size	Code	Compression	Colour
S	68310	15-20mmHg	Black	S	68315	20-30mmHg	Black
M	68311	15-20mmHg	Black	M	68316	20-30mmHg	Black
L	68312	15-20mmHg	Black	L	68317	20-30mmHg	Black
XL	68341	15-20mmHg	Black	XL	68318	20-30mmHg	Black
S	68363	15-20mmHg	Khaki	S	68349	20-30mmHg	Khaki
M	68326	15-20mmHg	Khaki	M	68321	20-30mmHg	Khaki
L	68327	15-20mmHg	Khaki	L	68322	20-30mmHg	Khaki
XL	68362	15-20mmHg	Khaki	XL	68323	20-30mmHg	Khaki

LADIES OPAQUE KNEE HIGH SOCK

Helps to energise tired legs, ankles, & feet. Improve circulation and prevent swelling. Featuring soft comfort non-restrictive bands whilst maintaining placement throughout the day. Reinforced heel/toe for durability, sheer two way stretch nylon & spandex.

Circular knitted, available in Ready-To-Wear sizes Sml - XX/Large, Ccl.I or Ccl.II, in either Sand/Black as standard, with other colours available upon special request.

SIZE CHART

Size	Ankle	Calf	Length
S	6.25"-7.25" (16-18cm)	11"-13.50" (28-34cm)	12"-18" (31-46cm)
M	7.50"-9" (19-22cm)	13"-15.50" (33-39cm)	13"-19" (33-48cm)
L	9.25"-10.75" (24-27cm)	15"-17.50" (38-44cm)	14"-20" (36-51cm)
XL	11"-11.50" (28-29cm)	17.75"-18.50" (45-47cm)	15"-22" (38-56cm)
XXL	11.75"-12.75" (30-32cm)	18"-19.50" (46-50cm)	15"-22" (38-56cm)

Size	Code	Compression	Colour	Size	Code	Compression	Colour
S	68128	15-20mmHg	Black	S	77108	20-30mmHg	Black
M	68124	15-20mmHg	Black	M	77104	20-30mmHg	Black
L	68125	15-20mmHg	Black	L	77105	20-30mmHg	Black
XL	68132	15-20mmHg	Black	XL	77106	20-30mmHg	Black
XXL	68121	15-20mmHg	Black	XXL	77107	20-30mmHg	Black
S	68200	15-20mmHg	Sand	S	75112	20-30mmHg	Sand
M	68201	15-20mmHg	Sand	M	75113	20-30mmHg	Sand
L	68202	15-20mmHg	Sand	L	75114	20-30mmHg	Sand
XL	68203	15-20mmHg	Sand	XL	75115	20-30mmHg	Sand
XXL	68204	15-20mmHg	Sand	XXL	75116	20-30mmHg	Sand

SMARTKNIT® SEAMLESS


by THERAFIRM®



Smartknit® seamless socks feature a non-binding Halo Top® that fits comfortably whilst minimising any tourniquet effect



Produced using high-stretch core spun fiber maximising products stretch whilst maintaining correct fit.



Moisture wicking CoolMax®* fibers provide dry comfortable coolness.
*Coolmax is a trademark of Invista



Black and Grey available with X-Static Silver Fiber Yarns: Anti-Bacterial, Odor Reducing, Temperature Central Fibers



Smooth seam free toe, with no pressure points. Ideal for diabetic patients



WHERE TO MEASURE

1. Measure the largest portion of the calf
2. Measure just above the ankle bone

KNIT-RITE DIABETIC SMARTKNIT SEAMLESS

DIABETIC SMARTKNIT SEAMLESS

SmartKnit® Seamless Socks by Therafirm® are produced using state-of-art machinery to totally eliminate seams thus risks associated with them. SmartKnit® Socks feel like a second skin and provide enhanced comfort and relief for those who have diabetes, arthritis, or sensitive feet.

SmartKnit® Seamless Socks are made from superior moisture wicking yarns including X-Static SILVER fibers to wick moisture away from the skin to the exterior of the sock, moving moisture away from the foot. These technical yarns result in a sock which effectively controls moisture and helps reduce blister-causing friction.

SmartKnit® Seamless Socks patented knitting process use hi-tech, high-stretch Core-Spun yarns to ensure a sock which 'hug' you feet to ensure maximum protection to sensitive feet, whilst not being bulky inside a shoe when worn. The non-binding Halo Top feature ensures that whilst the sock stays where required they minimise any grip or indentations and offer a very comfortable fit.

All Black or Grey SmartKnit® Seamless Socks contain industry leading X-Static Yarns. This use of high-tech Silver Fibers naturally reduce foot odor, wick away perspiration to keep the foot dry, and maintain the foot temperature by conducting both heat and cold.

SIZE CHART

Size	Men's shoe size	Womans
S	5-7	5-8
M	8-10	9-11
L	11-13	12-14
XL	14+	-

CREW LENGTH		
Size	Code	Colour
S	71122	White
M	71132	White
L	71142	White
XL	71152	White
S	71722	Black
M	71732	Black

WIDE CREW LENGTH		
Size	Code	Colour
S	71162	White
M	71172	White
L	71182	White
XL	71192	White

OVER CALF (BK)		
Size	Code	Colour
S	73122	White
M	73132	White
L	73142	White
XL	73152	White
S	73722	Black
M	73732	Black
L	73742	Black
XL	73752	Black


Size	Code	Colour
L	71742	Black
XL	71752	Black
S	71622	Grey
M	71632	Grey
L	71642	Grey
XL	71652	Grey

Size	Code	Colour
S	71762	Black
M	71772	Black
L	71782	Black
XL	71792	Black

Size	Code	Colour
S	73422	Navy
M	73432	Navy
L	73442	Navy
XL	73452	Navy



LYMPHAPRESS PUMP ADVANTAGES



24 cell operations (12 per limb).



Single limb operation (12 cells) using supplied plug.



Single quick connect connectors for hose bundles.



Pretherapy modes (optimal only).



Delivers fully variable 20-90mmHg pressure.



User friendly programming & quiet operation.

UNIQUE GARMENTS




4, 8, 12 or 24 Overlapping cells, garment dependent. Ensuring consistent pressure gradient along limb.



Manufactured using state of the art High frequency welding technology.



Integrated hose bundles. No need to connect hoses.



Available in a range of styles and sizes to ensure adequate coverage and fit.

PUMPS & EQUIPMENT

(4/8 Cell)

LYMPHA PRESS® MICRO (STD)



Features

- LymphaPress® Micro is a sequential compression device designed for use by either the patient at home or by practitioners in clinics.
- Automatically assesses and adjusts the garment fill and pause interval time according to the garment being used.
- Complete garment deflation each cycle for patient comfort and maximum vascular refill.
- Easy to use.
- Lightweight and quiet.
- Dual Voltage - Ideal for travel.



Small compact design.



Easily lockable compression setting dial.



Easy to use timer and compression settings dials.



4 chamber or 8 chamber garments available.



Single or dual garment outlet.



Lightweight and quiet operation



Dual Voltage - Ideal for use overseas thus continuing treatment.

Technical Specifications

Cycle	Sequential
Pressure Range	20-80 mmHg
Size	12 x 22 x 17cm
Weight	1.5kg
Electrical	230V / 115V 50/60 Hz
Code	L10015G4

LYMPHA PRESS® MICRO (OPTIMAL)



Based on the reliable and durable Lympha Press® Micro line, the Lympha Press Micro (Optimal) is a pressure control device with the added features of both Pre and Wave Therapy treatments.

For use within both the clinical and home environment. The Pretherapy Treatment is a gentle massaging of the proximal lymphatics, with the purpose of easing the lymph liquids transit, based on MLD (manual lymph drainage) procedures. The Pretherapy starts the treatment at the proximal end of the treated limb. Lympha Press® Micro (Optimal) then moves on to the main treatment, either sequential buildup pressure or Wave Peristaltic massage.

Indications

- Lymphoedema of legs and arms.
- Post-mastectomy arm Lymphoedema.
- Healing of venous stasis ulcers.
- Promotion of wound recovery.
- Relief from muscular and nervous pain.
- Relief of discomfort from leg fatigue.

Additional Features to the Lympha Press Micro (Optimal)

- Fully variable pressure range 20 to 80mmHg.
- Sequential Compression and Peristaltic Compression Therapy Cycles.
- Pretherapy and Post Therapy available with all Compression Cycles.
- Variable Pre Therapy Region: Distal, Medial and Proximal.
- Variable Post Therapy Region: Distal, Medial and Proximal.
- Adjustable Pause time between compression cycles.
- Adjustable hold time at end of compression cycles.

Technical Specifications

Cycle	Sequential Multi Cycles
Pressure Range	20-80 mmHg
Size	12 x 22 x 17cm
Weight	1.5kg
Electrical	230V / 115V 50/60 Hz
Code	L10015H1-AU

PUMPS & EQUIPMENT

(12 Cell)

LYMPHA PRESS® MINI PRESS



The new Lympa Press® Mini Press, model 960 offers fully pressure controlled 12 chamber sequential compression, with user selected treatment time and pressure from 20-80 mmHg. The product is easy to operate, Multi-Voltage and ideal for treatment of a wide variety of patients at home or in the clinic.

Features

- Sturdy and reliable , user friendly for clinic and home use
- 12 chambers cost effective device
- Pressure Controlled
- Multi Voltage
- Compatible with Lympa Press ® famous range of unique and innovative 12 chamber garments

Indications

Primary Lymphoedema

- Congenital (Milroy's disease)
- Proximal obliteration
- Distal obliteration

Secondary Lymphoedema

- Post mastectomy
- Groin dissection
- Chronic edema of venous origin
- Post traumatic edema
- Lymphatic filariasis

Other Indications

- Lipedema
- Venous insufficiency
- Varicose veins
- Leg ulcers
- Promotion of wound recovery
- Postoperative venous ligation or stripping

Contraindications

Use of the Mini Press® system is not recommended in the presence of one or more of the following conditions:

- Known or suspected deep vein thrombosis (DVT) or pulmonary embolism
- During the inflammatory phlebitis process
- Acute infection of the affected limb
- Decompensated cardiac failure
- Severe arteriosclerosis or other ischemic vascular disease
- Any circumstance where increased venous and lymphatic

Due to movement of fluids in the body when using the device, exercise caution when using on patients with heart disease. High pressure is not recommended for patients who have peripheral occlusion disease. The abdominal area should not be treated during pregnancy.

Technical Specifications

Cycle	Sequential
Pressure Range	20-80 mmHg
Size	34cm x 24cm x 19cm
Weight	2.9 kg
Electrical	100-240 VAC 50-60Hz
Code	L10015H1-AU

LYMPHA PRESS® OPTIMAL PLUS



Lympha Press® Optimal Plus, model 912 is the most versatile and advanced dynamic compression system for lymphoedema treatment on the market today.

It offers unique, state of the art features and treatment modes including:

- Lympha sequential mode for effective reduction of edema and fibrosis.
- Wave peristaltic mode for comfortable, effective treatment in patients with distal discomfort.
- Pretherapy™ designed on the principles of manual lymph drainage
- Post-therapy™ additional therapy to a specific area after the main treatment is completed.

Features

- Adjustable range 20 to 90 mmHg in a gradient, to enhance distal to proximal flow.
- Automatic adjustable Pause time to ensure vein refill
- User-friendly, clear interface
- Programmable via computer, for pressure adjustment of individual chambers
- Compatible with comprehensive range of unique and innovative garments
- Air deflation after treatment
- Full compliance meter. Easy data download via USB port.
- Multi Voltage
- Cost Effective

Indications

The device is intended for use by medical professionals, and patients who are under medical supervision, in treating many conditions such as:

- Primary lymphoedema
- Secondary lymphoedema
- Lipedema
- Venous insufficiency
- Venous stasis ulcers
- Dysfunction of the muscle pump
- Post mastectomy edema
- Edema following trauma and sports issues
- Reduction of wound healing time
- Reduction of pain and swelling after injury and surgery

Contraindications

Use of the Lympha Press Optimal® Plus system is not recommended in the presence of one or more of the following conditions:

- Known or suspected deep vein thrombosis (DVT) or pulmonary embolism
- During the inflammatory phlebitis process
- Acute infection of the affected limb
- Decompensated cardiac failure
- Severe arteriosclerosis or other ischemic vascular disease
- Any circumstance where increased venous and lymphatic return is undesirable

Due to movement of fluids in the body when using the device, exercise caution when using on patients with heart disease. High pressure is not recommended for patients who have peripheral occlusion disease. The abdominal area should not be treated during pregnancy.

Technical Specifications

Cycle	Sequential, Wave, Pretherapy, Post-therapy
Pressure Range	20-90 mmHg
Size	34cm x 24cm x 19.8cm
Weight	4.5 kg
Electrical	100-240 VAC 50-60Hz
Code	L10008C1 - AU

*Also available with Bluetooth connection

Code	L10008C2 - AU
------	---------------

PUMP GARMENTS

8 CELL COMFY ARM SLEEVE WITH INTEGRATED HOSE



The 8 Cell Comfy Arm Sleeve offers an ideal way to treat Lymphoedema in the Upper limbs.


The treated area does not stop at the upper arm but also includes the upper shoulder, back, chest and abdomen.

For use with Lympha Press Micro (Std) and Lympha Press Micro (Optimal).

Code	Description	Torso Dimensions
L41LT100	Comfylite Arm 8	Available in one size. Sizing is dependent on fastener type used
L42LOFSL	Left Arm Fastener - Small	
L42LOFML	Left Arm Fastener - Medium	
L42LOFLL	Left Arm Fastener - Large	
L42LOFSR	Right Arm Fastener - Small	
L42LOFMR	Right Arm Fastener - Medium	
L42LOFLR	Right Arm Fastener - Large	



8 Cell Comfy Arm Sleeve with Integrated Hose Sizing Chart



Notes: All dimensions are in centimeters (cm) General tolerance: +2 cm

12 CELL BI-LATERAL LYMPHA® JACKET



Designed to give best lymph drainage for both arms and upper abdomen.


Very useful and efficient where both or only one arm suffers from Lymphoedema.

For use with Lympa Press Mini and Lympa Press Optimal Plus.

Code	Style	Arm 1	Arm 2	Chest Circ.
L4MJ2125	Short	800mm	570mm	1050-1230mm
L4MJ2135	Long	960mm	720mm	1050-1230mm



Lympa Sizing Chart



12 CELL COMFY® ARM SLEEVE



The Comfy Arm Sleeve 1 - 75 offers a new and greatly improved way to treat Lymphoedema in the upper limbs.


The treated area does not stop at the upper arm but now includes the armpit as well as the upper shoulder, one side of the back, chest and abdomen.

For use with Lympa Press Mini Press and Lympa Press Optimal Plus.

Code	Description	Torso Dimensions
L41L0175R	Comfy Arm Sleeve, Right Arm and Torso	Available in one size. Sizing is dependent on fastener type used
L41L0175L	Comfy Arm Sleeve, Left Arm and Torso	Available in one size. Sizing is dependent on fastener type used
L42LOFSL	Left Arm Fastener - Small	
L42LOFML	Left Arm Fastener - Medium	
L42LOFLL	Left Arm Fastener - Large	
L42LOFSR	Right Arm Fastener - Small	
L42LOFMR	Right Arm Fastener - Medium	
L42LOFLR	Right Arm Fastener - Large	



12 Cell Comfy Arm Sleeve with Integrated Hose Sizing Chart



PUMP GARMENTS

4 CELL ARM SLEEVE



The 4 Cell Arm Sleeve offers patients an ideal and easy to don option for treating Upper Limb Lymphoedema.

The 4 Cell arm sleeve has an axilla cut out so that the shoulder is treated along with the arm.


For use with Lympha Press Micro (Std) and Lympha Press Micro (Optimal).

Code	Description
L41F1290H	4 Cell Arm Sleeve



4 Cell Arm Sleeve Sizing Chart

Part No.	Model/Size	A (Upper Arm)	C (Distal End)	Length
L41F1290H	LP4ARM M/L	58	41	88
L415A087	Expander C	Adds 10 cm to the circumference of the sleeve		



12 CELL ARM SLEEVE



Special design makes arm and shoulder treatment comfortable and easy.

The new 12 Cell arm-shoulder sleeve provides therapy to the entire shoulder area as well as the arm.


The design allows treatment for either right or left arm with a single garment.

For use with Lympha Press Mini Press and Lymph Press Optimal Plus.

Code	Description
L4111C2	12 Cell Arm Sleeve - Med
L4111C3	12 Cell Arm Sleeve - Large



12 Cell Arm Sleeve Sizing Chart



This garment has two zippers, for adjustment of the circumferential fit. The measurements indicated are the maximum (outmost zipper)

Part No.	Model/Size	Circumference					Length		
		A (Upper Arm)	B	C (Distal End)	Shoulder	Strap (S)	L1	L2	1/2 L2
L4111C2	LP12 ARM M	63	49	41	76	92 - 156	91	72	36
L4111C3	LP12 ARM L	74	56	49	88	92 - 165	106	91	40

PUMP GARMENTS

4 CELL BARIATRIC LEG SLEEVE



Bariatric Leg Sleeve is a compression garment for the treatment of full legs - foot, calf and thigh.

For use with Lympha Press Micro (Std) and Lympha Press Micro (Optimal).

Code

L41F2365H

L41F2375H

L41F2385H

8 CELL COMFY® LEG SLEEVE



Comfy Leg Sleeve is a compression garment for the treatment of full legs - foot, calf and thigh.

For use with Lympha Press Micro (Std) and Lympha Press Micro (Optimal)

Code


L41H2365

L41H2375

L41H2385



Leg Sleeve Sizing Chart



Garment P/N	Description	Dimensions								Type	Expander P/N	Dimensions with Expander T						
		A	B	C	D	E	F	L	A			B	C	D	E	F	L	
		Thigh	Above Knee	Calf	Ankle	Instep	Foot	Inseam			Thigh	Above Knee	Calf	Ankle	Instep	Foot	Inseam	
L41F2465H	Lympha Press 4 4-65	105	98	86	67	60	41	65	FT	L415T465	119	112	99	72	64	41	65	
L41F2475H	Lympha Press 4 4-75	107	99	86	72	60	41	75	FT	L415T465	121	113	99	77	64	41	75	
L41F2485H	Lympha Press 4 4-85	109	105	93	67	60	43	85	GT	L415T485	123	118	106	71	64	43	85	
L41H2365	Lympha Press 8 3-65	78	71	65	50	54	38	65										
L41H2375	Lympha Press 8 3-75	81	74	66	50	50	38	75										
L41H2385	Lympha Press 8 3-85	86	80	68	53	51	38	85										

Notes: All dimensions are in centimeters (cm) General tolerance: +2 cm

12 CELL COMFY® LEG SLEEVE




Comfy Leg Sleeve is a multiple air-cell compression garment for the treatment of full legs - foot, calf & thigh.

For use with Lympha Press Optimal Plus and Lympha Press Mini Press.

Code	Size	A	B	C	D	E	F	L
L41L2285	2-850mm	700mm	600mm	470mm	420mm	480mm	260mm	840mm
L41L2295	2-950mm	780mm	680mm	510mm	430mm	480mm	260mm	950mm
L41L2365	3-650mm	810mm	730mm	600mm	560mm	540mm	370mm	650mm
L41L2375	3-750mm	850mm	760mm	610mm	560mm	540mm	370mm	750mm
L41L2385	3-850mm	900mm	800mm	630mm	560mm	540mm	370mm	850mm

1. General tolerance: +2 millimeters. 2. A, B, C, D & E dimensions are circumferences.

COMFY® LEG SLEEVE PARALLEL EXPANDER



Comfy Leg Sleeve Parallel Expander allows an increase in the circumference of the Comfy Leg Sleeve.

Each expander adds 100mm from the top to the bottom of the sleeve.


Code	Use with Sleeve
L415A097	L41L2385
L415A087	L41L2285, L41L2375
L415A101	L41L2295
L415A078	L41L2265

Comfy Leg Sleeve Details		Sleeve Dimensions							Parallel Expander						
Cat No.	Size	A	B	C	D	E	F	L	Cat No.	Size	Sleeve Dimensions with Expander				
											A	B	C	D	E
L41L2285	2-85	70	60	47	43	48	26	84	L415A087	C	80	70	57	53	58
L41L2295	2-95	80	72	60	43	48	26	95	L415A101	E	90	82	72	58	58
L41L2365	3-65	81	73	60	56	54	37	65	L415A078	B	91	83	70	66	64
L41L2375	3-75	85	76	61	56	54	37	75	L415A087	C	95	86	71	66	64
L41L2385	3-85	90	80	64	56	54	37	85	L415A097	D	100	90	74	66	64
L41L2475	4-75	76	-	-	40	-	-	87	L415A087	C	86	-	-	50	-

1. All dimensions are in centimeters. 2. General tolerance: ± 2 centimeters. 3. A, B, C, D, E dimensions are circumferences.

PUMP GARMENTS

COMFY® LEG SLEEVE TAPERED EXPANDER



Comfy Leg Sleeve Tapered Expander allows an increase in the circumference of the Comfy Leg Sleeve.

Each expander adds 100mm above the knee tapered to 40mm at the ankle.

Code	Use with Sleeve
L415T097	L41L2385
L415T087	L41L2285, L41L2375
L415T101	L41L2295
L415T078	L41L2265

Comfy Leg Sleeve Details		Sleeve Dimensions							Parallel Expander						
Cat No.	Size	A	B	C	D	E	F	L	Cat No.	Size	Sleeve Dimensions with Expander				
											A	B	C	D	E
L41L2285	2-85	70	60	47	43	48	26	84	L415T087	CT	80	70	56	49	52
L41L2295	2-95	80	72	60	43	48	26	95	L415T101	ET	90	82	60	49	52
L41L2365	3-65	81	73	60	56	54	37	65	L415T078	BT	91	83	70	59	60
L41L2375	3-75	85	76	61	56	54	37	75	L415T087	CT	95	86	71	59	60
L41L2385	3-85	90	80	64	56	54	37	85	L415T097	DT	100	90	73	59	60
L41L2475	4-75	76	-	-	40	-	-	87	L415T087	CT	86	-	-	44	-

1. All dimensions are in centimeters. 2. General tolerance: ± 2 centimeters. 3. A, B, C, D, E, F, L dimensions are circumferences.

LYMPHA® PANTS EXPANDER - FOR USE ON ALL BI-LATERAL PANTS



Lympha Pants Expanders® allow extending the circumference of Lympha Pants®.

Use Lympha Pants Expanders® with Mego Afek: Lympha Pants®. You can use up to three pants expanders on each pant leg. Each expander (right & left) adds about 10cm to the circumference of the pant (leg & abdomen). Pants Expanders come in sets of two. One for each side of the pants.

Instructions for use

1. Open the zipper on the pant leg that you want to expand.
2. Lay out the Pant Expander® next to one side of the open zipper.
3. Close one side of the Pant Expander® zipper with one side of the pants zipper.
4. Close 1/3 of the other side of the Pant Expander® zipper with the other side of the pants zipper.
5. Use the expanded Lympha Pants® in the standard manner.

Code	Unit
L4Z5LP00	Pair

12 CELL BI-LATERAL LYMPHA PANTS



Bi-Lateral Lympa Pants are designed to offer optimal lymph drainage for both legs and the lower abdomen.

Bi-Lateral Lympa Pants are also very effective in cases where only one leg suffers from Lymphoedema.

Use Bi-Lateral Lympa Pants with the following Mego Afek Devices: Lympa Press® Optimal & Lympa Press® Mini

Instructions for wearing the Bi-Lateral Lympa Pants

1. Arrange the Bi-Lateral Lympa Pants on the surface where treatment will take place.
2. Open the main zipper 2/3 & lie down inside the Bi-Lateral Lympa Pants.
3. There are three zippers for three different pant sizes. Close the zipper which adapts the size to meet your needs, closing the pants.
4. If the pants are still too small, add Expanders® to adapt the pants size to your size.
5. On each leg, near the ankle, there are three (3) small zippers. Open the one that is closest to the patient's heel. It will allow the foot to remain as close as possible to its normal position (90° to leg) during treatment.

Recommendations

- Use Bi-Lateral Lympa Pants while lying down
- Wear cotton training pants under the Bi-Lateral Lympa Pants

Connecting the hose bundle

Insert the hose bundle connectors (in direction pointed by red arrows) into the corresponding air socket receptacles of the device.


Disconnecting the hose bundle

Press down the snap tooth to release the hose bundle and gently pull out the bundle from the device.

Code	Size
L4MP2136C	Regular
L4MP2226C	X-Large



Bi-Lateral Lympa Pants Sizing



PUMP GARMENTS

8 CELL BI-LATERAL PCD LYMPHA PANTS



Bi-Lateral PCD Lympa Pants are designed to offer optimal lymph drainage for both legs and the lower abdomen.

Bi-Lateral PCD Lympa Pants are also very effective in cases where only one leg suffers from Lymphoedema.

Use Bi-Lateral PCD Lympa Pants with the following Mego Afek Devices:
Lympa Press® Micro (Optimal)

Instructions for wearing the Bi-Lateral PCD Lympa Pants

1. Arrange the Bi-Lateral PCD Lympa Pants on the surface where treatment will take place.
2. Open the main zipper 2/3 & lie down inside the Bi-Lateral PCD Lympa Pants.
3. There are three zippers for three different pant sizes. Close the zipper which adapts the size to meet your needs, closing the pants.
4. If the pants are still too small, add Expanders® to adapt the pants size to your size.
5. On each leg, near the ankle, there are three (3) small zippers. Open the one that is closest to the patient's heel. It will allow the foot to remain as close as possible to its normal position (90° to leg) during treatment.

Recommendations

- Use Bi-Lateral PCD Lympa Pants while lying down
- Wear cotton training pants under the Bi-Lateral PCD Lympa Pants

Connecting the hose bundle

Insert the hose bundle connectors (in direction pointed by red arrows) into the corresponding air socket receptacles of the device.

Disconnecting the hose bundle


Press down the snap tooth to release the hose bundle and gently pull out the bundle from the device.

Code

L4MPH004



Bi-Lateral PCD Lympa Pants Sizing



8 CELL BI-LATERAL PCD LITE LYMPHA PANTS



Bi-Lateral PCD Lite Lympa Pants are designed to offer optimal lymph drainage for both legs and the lower abdomen.

Use Bi-Lateral PCD Lite Lympa Pants with the following Mego Afek Devices: Lympa Press® Micro (Optimal)

Instructions for wearing the Bi-Lateral PCD Lite Lympa Pants

1. Arrange the Lympa Pants (Leg Specific) on the surface where treatment will take place.
2. Open the main zipper 2/3 & lie down inside the PCD Lite Lympa Pants.
3. If the pants are still too small, add Expanders® to adapt the pants size to your size.
4. On each leg, near the ankle, there are three (3) small zippers. Open the one that is closest to the patient's heel. It will allow the foot to remain as close as possible to its normal position (90° to leg) during treatment.

Recommendations

- Use Bi-Lateral PCD Lympa Pants while lying down
- Wear cotton training pants under the Bi-Lateral PCD Lite Lympa Pants

Connecting the hose bundle

Insert the hose bundle connectors (in direction pointed by red arrows) into the corresponding air socket receptacles of the device.

Disconnecting the hose bundle

Press down the snap tooth to release the hose bundle and gently pull out the bundle from the device.


Code

L4MPH003 - Regular

L4MPH005 - Small



Bi-Lateral PCD Lite Lympa Pants Sizing



PUMP GARMENTS

LYMPHAPOD®



LymphaPod® is the only pneumatic compression system designed specifically for the morbidly obese.

- Revolutionary, innovative construction.
- Designed especially for patients with limited mobility.
- Unique “sleeping bag” design for easy donning and doffing.
- One size fits all, from 150 to 300+ kilograms.

Code

L4MP800L



TAPERED LYMPHAPOD® EXPANDER



LymphaPod® allows extending the circumference of the LymphaPod®.

Each expander adds 170mm at the abdomen tapered to 90mm at the foot.

Code	Side
L41580AL	Left
L41580AR	Right

MEGO AFEK: LYMPHA PRESS

PROTOCOLS FOR LYMPHA PRESS® THERAPY

Treatment for lymphoedema, venous insufficiency, venous stasis ulcers, and dysfunction of the “muscle pump.”

Australian TGA restricts this device to sale and use by, or on the order of, a licensed medical practitioner. Without referrals devices contained in this catalogue cannot be purchased

Special Considerations: Congestive heart failure

Decompensated cardiac failure (CHF) is a contraindication to pump use. Patients who have had recent hospitalization for CHF should receive physician clearance prior to pump use. Patients with a history of CHF should be monitored for symptoms of shortness of breath, increased rate or labored breathing, lung crackles, cardiac gallop, and jugular vein distension. Patients with CHF who avoid lying flat due to shortness of breath should have close physician monitoring if they use pneumatic compression.

Monitor symptoms during and after treatment. Start treatment at low pressure (<35 mmHg) and short (<30 minute) treatment and gradually increase pressure and time according to tolerance and treatment response. For initial leg treatments in patients with CHF, start by using pump on one leg, and then follow up with the second leg the next day or treatment session. If after pumping there is no worsening in symptoms, simultaneous bilateral treatment can be applied. If pants-type appliance garments are used, treat bilaterally at low pressures for short periods of time (30 minutes, <35 mmHg) and increase time and pressure gradually per treatment session, according to tolerance.

Patients undergoing treatment in hospital can receive IV diuretics prior to treatment. Start treatment after urine output is established. Pulse oximetry can be used to monitor lung function.

Mixed venous and arterial disease

Some patients with venous stasis ulcers may have arterial insufficiency as well as venous insufficiency. Use Lympha Press® and The Petite Basic System™ with caution (lower pressures, frequent inspection of skin) on patients with ABI from 0.6 to 0.8. Lympha Press® and The Petite Basic System™ pumps should not be used on patients with ABI <0.6.



Deep vein thrombosis (DVT)

Patients with a history of DVT should have a Doppler ultrasound to rule out DVT prior to initiating treatment.

Cellulitis infection

Patients should avoid pneumatic compression treatment during an acute infection; treatment may be resumed 48-72 hours after beginning antibiotics.

Cellulitis infection

Patients who are unable to operate pump, don, doff, or adjust the appliance garments should have an appropriate family member or other caregiver available to assist throughout the treatment session.

QUICK REFERENCE GUIDE

This document offers guidance for the licensed medical practitioner in selection of treatment parameters. These treatment parameters (pressure, frequency, and duration of treatment) should be selected according to individual patient characteristics, including edema, fibrosis, wound status, body mass and response to a trial treatment session.

Patient Evaluation

Obtain medical history: Surgeries, co-morbidities; pain/ sensation levels, functional ability/mobility, edema, and skin condition including fibrosis, scars and wounds. Obtain baseline (pre-treatment) oedema measurements.

Rule out general contraindications to pump use: Known or suspected deep vein thrombosis or pulmonary embolism; inflammatory phlebitis process; acute infection of the affected limb; decompensated cardiac failure; severe arteriosclerosis or other ischemic vascular disease; any circumstance where increased venous and lymphatic return is undesirable. Due to movement of fluids in the body when using the device, be careful with patients who have heart disease. High pressure is not recommended for patients with peripheral vascular disease. The abdominal area should not be treated during pregnancy.

Perform a trial to determine response to treatment time and pressure

Follow the guidelines for pressure selection provided in Tables 1-3. Begin with the lower pressure ranges, and adjust pressure upward within the recommended ranges according to patient tolerance and response. Increase pressure if no adequate response is seen. The patient should feel comfortable during treatment. If patients are insensate or have fragile skin, begin at conservative levels of compression and observe for skin changes. If using an appliance garment that covers the abdomen, set the pressure at this level below 40 mmHg.

If patient feels discomfort, investigate source. Reducing the compression level may alleviate discomfort. If discomfort persists, stop treatment and determine whether pump treatment should be continued.

Stop treatment and assess patient if the patient experiences dizziness, numbness, tingling, cramping, coldness, extremity or chest pain, or shortness of breath during or immediately following treatment. After assessment, if treatment is continued, compression levels and/or treatment time should be resumed at a reduced level.

After completion of the trial treatment session, check for edema reduction, softening and tissue extensibility changes. Inspect the skin to ensure there are no adverse effects.

Prior to use

1. Rule out contraindications.
2. Trial pump on patient, starting with lower pressure and increasing pressure if no adequate response seen. Patient should feel comfortable during treatment. If patient feels discomfort, investigate source. Reducing the compression level may alleviate discomfort.

Table 1: Venous Ulcer/Venous Stasis			
Compression Level	Pressure	Appropriate for:	Frequency/ Duration
Low compression	<40 mmHg	Patients with wounds, insensate, or fragile tissue (can increase with monitoring if no adverse response seen)	60 minutes once or twice per day. Shorter, more frequent sessions may be helpful for patients who can comply. Evening or afternoon sessions will provide the most response.
Moderate compression	40-60 mmHg	Most patients	
High compression	>60mmHg	Extreme fibrotic changes; lipodermatosclerosis	
Table 2: Lower Extremity Lymphoedema			
Compression Level	Pressure	Appropriate for:	Frequency/ Duration
Low to moderate compression	30-60 mmHg	Low to normal body mass, Stage 1 or Stage 2 lymphoedema, patients with wounds, insensate or fragile tissue (can increase with monitoring if no adverse response seen)	60 minutes once or twice per day. Shorter, more frequent sessions may be helpful for patients who can comply. Evening or afternoon sessions will provide the most response.
Moderate to high compression	>60 mmHg	Stage 3 lymphoedema, heavy fibrosis and/or high body mass. Use <40 mmHg compression over abdomen for comfort	
Table 3: Upper Extremity Lymphoedema / Post-Mastectomy Lymphoedema			
Compression Level	Pressure	Appropriate for:	Frequency/ Duration
Low to moderate compression	30-50 mmHg	Low to normal body mass, Stage 1 or Stage 2 lymphoedema, patients with wounds, insensate or fragile tissue (can increase with monitoring if no adverse response seen)	60 minutes once or twice per day.
Moderate to high compression	40-60 mmHg	Stage 2 or Stage 3 lymphoedema, heavy fibrosis, or high body mass. Use <40 mmHg compression over torso for comfort	

After use

After the pumping session, inspect for oedema reduction, softening and tissue extensibility changes. Assess the skin for any adverse effects (skin breakdown, irritation, or other).



ONE DEVICE MULTIPLE USES

Sports Therapy

Versatile solutions for a range of sports injuries, with working on site possibilities. Aiding in maintenance and recovery.

Oncology and Lymphoedema

3x more efficient than Manual Lymphatic Treatment.

Orthopaedics

Reduce (Pre and Post Op) swelling for earlier patient discharge.

MSD

Localised and comprehensive treatment of musculoskeletal conditions.

Neurology Pain Management

Non-pharmacologic approach to treatment by concentrating on fascial structures and trigger points.



Watch Lympha touch in action



Ergonomic handle for patient and therapist comfort




Battery or main use



Sports mode available for up to 350mm Hg pressure



Adjustable pulsation and pressure control




In-built treatment protocols




Lightweight and portable for on site use

OPERATING PRINCIPLE



Tissue Expansion

Negative pressure expands the tissue, creating movement of the anchoring filaments and the endothelial holes expand.



Endothelial expansion

Metabolic waste moves into the lymphatic capillary network - lymph is transported through these capillaries to the pre-collectors.



Lymph Flow Increased

Negative pressure further enhances the flow of lymph into the collectors (the body's natural pumping mechanism).

Code	Description
LT-1001	LymphaTouch, Negative Therapy Pressure Device
LT-1200	LymphaTouch, Sports Upgrade - Negative Pressure Up to 350mmHg
LT-1901	LymphaTouch Stand
LT-1004-10	Treatment Cup, 10mm, 4 Piece
LT-1004-35	Treatment Cup, 35mm, 4 Piece
LT-1004-50	Treatment Cup, 50mm, 4 Piece
LT-1004-60	Treatment Cup, 60mm, 4 Piece
LT-1004-80	Treatment Cup, 80mm, 4 Piece
LT-1004-MIX	Treatment Cup, 1x35mm, 1x50mm, 1x60mm, 1x80mm
LT-1315-5	LymphaTouch Filters for LT-1001, 500pcs

Absolute Contraindications

- Acute Deep Vein Thrombosis
- Acute Infection
- Congestive Heart Failure, Cardiac Oedema
- Kidney Dysfunction
- Conditions where increased venous or lymphatic return is undesirable

Relative Contraindications

- Active Cancer - seek medical advice in advance
- Pregnancy - seek medical advice in advance
- Addition: All contraindications for MLD/Physiotherapy Apply
- Severe ankle and leg swelling
- Sensitivity disorders of the extremities



80mm



60mm



50mm



35mm



10mm



an **OPED** company

OPC Health
26-32 Clayton Rd, Clayton VIC 3168
T: 1300 672 937 (1300 OPC YES)
E: sales@opchealth.com.au
Visit our website at
opchealth.com.au

*OPC Health reserves the right to make any amendments and/or variations to product or pricing without notice.