General Information	Ship To:	Special Instructions:
Ordered by :	Address:	
•	City:	
Phone:	State: Zip Code: Country:	
BILLING: P.O. Number	Phone:Attention:	
Townsend Account Number:	FIIONEAttention.	
Bill To:	Fit Date: If known, please indicate the date you are	
Address:	scheduled to fit the patient:	
City:	Shipping Preference:	
State: Zip Code: Country:	Ground 2-Day P.M. 2-Day A.M.	
	Next Day P.M. Next Day A.M.	
Phone: Fax:	(If no preference is indicated, this order will be shipped 2 Day P.M.)	
	Note: We do not ship new or repaired products directly to patients.	
FAB SERIES This is a unique elastic polyester blend casting tape (not fiber glass, not plaster) that is easier to contour than traditional fiber glass. It has the molding characteristics of plaster without the mess. FAB2 (2" X 4 yards; white tape) Boxes:	FG SERIES This is a multi-purpose fiber glass casting tape for use/wear by a patient for stabilization or positional support following an injury or surgery. The FG series is also suitable for creating patient molds for use in the fabrication of orthopedic braces and prosthetic devices. FG2 (2" X 4 yards; cast tape) Color: White Red Lime Green Boxes: 1 2 3 4 Other FG3 (3" X 4 yards; cast tape) Color: White Red Lime Green Boxes: 1 2 3 4 Other FG4 (4" X 4 yards; cast tape) Color: White Red Lime Green Boxes: 1 2 3 4 Other FG4 (4" X 4 yards; cast tape) Color: White Red Lime Green Boxes: 1 2 3 4 Other	Instructions For Use NOTE: KEEP TAPE WET Tear Top Off Of Bag, Leave Roll In Bag, & Immerse In Warm Water Squeeze Out Excess Water (With Roll In Bag), Submerge/ Squeeze Again Remove Roll From Bag, Roll Onto Limb & Contour Allow To Dry Until Rigid & Use Shark-
		Rigid & Use Shark- Nose Scissors
FAB5 (5" X 4 yards; white tape)	FG5 (5" X 4 yards; white tape)	To Remove
Boxes: 1 2 3 4 Other	Boxes:	

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313 Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722

TOWNSEND'S SHIPPING DEPARTMENT USE ONLY

RECEIVED ____

NEW BRACE SERVICE -- Original Brace Returned? Yes No

__ SHIPPED VIA _