

# ANTERIOR POSTERIOR TAPERED APDT-TP | APFR-TP

Gel Liner – With or Without  
Distal Attachment



## DESCRIPTION


The ALPS® AP Tapered Liner is gradually tapered from the anterior to the posterior to provide superior comfort. This liner is available in a pin and lock system or suction suspension. The AP Tapered Liner features our black high performance fabric with gel to assist in reducing bunching in the popliteal region.

PART #	SUSPENSION	L-CODE
APDT	LOCKING	L5673
APFR	CUSHION	L5679

Recommended for all activity levels.




## TOP FEATURES



- ✓ Constructed using high performance black fabric with gel to reduce bunching in the popliteal region
- ✓ Tapered – 6mm Anterior to 3mm Posterior


# CROSS SECTION DIAGRAM

6.0 MM TO 3.0 MM TAPERED





## HOW TO SIZE

Measurements taken 6cm from the distal end



## SIZING CHART

MEASURED SIZE (CM)	 LOCKING 6.0 MM or 3.0 MM TAPERED	 CUSHION 6.0 MM or 3.0 MM TAPERED
20-23 cm	APDT <b>20</b> - TP	APFR <b>20</b> - TP
24-25 cm	APDT <b>24</b> - TP	APFR <b>24</b> - TP
26-27 cm	APDT <b>26</b> - TP	APFR <b>26</b> - TP
28-31 cm	APDT <b>28</b> - TP	APFR <b>28</b> - TP
32-37 cm	APDT <b>32</b> - TP	APFR <b>32</b> - TP

\*Expired patent owned by ALPS and/or one of its subsidiaries



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862  
www.easyliner.com | info@easyliner.com  
2895 42nd Avenue North St. Petersburg, Florida 33714

